- ELEKTRODENPOSITIONERUNGSANLEITUNG
- ELECTRODE PLACEMENT GUIDE
- GUIA PARA LA COLOCACION DE ELECTRODOS
- GUIDE DE POSITIONNEMENT DES ÉLECTRODES
- GUIDA AL POSIZIONAMENTO DEGLI ELETTRODI
- OVERZICHT PLAATSING VAN DE ELEKTRODEN
- ELEKTRODEPLASSERINGSGUIDE
- ELEKTRODPLACERINGSGUIDE
- ELEKTRODIEN ASETTELU

CE 2797

MD

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15. REFERENZEN

Bundsen P, Peterson L-E, Selstam U

Pain relief in labor by transcutaneous electrical nerve stimulation. A prospective matched study.
Acta Obst Gyn Scand 1981; 60: 459-468

Grim L, Morey S

Transcutaneous electrical nerve stimulation for relief of parturition pain. A clinical report.
Physical Therapy 1985; 65: 337-340

Harrison R, Woods T, Shore M, Mathews G, Unwin A

Pain relief in labour using transcutaneous electrical nerve stimulation (TENS).
A TENS/TENS placebo controlled study in two parity groups.
Br J Obst and Gyn 1986; 93: 739-746

Kaplan B, Rabinerson D, Lurie S, Bar J, Krieser UR, Neri A

Transcutaneous electrical nerve stimulation (TENS) for adjuvant pain-relief during labor and delivery
Int Journal of Gynecology and Obstetrics 1998; 60: 251-255

Holgersson S, Johansson E

The effectiveness of TENS during labour – a survey
Research project in physiotherapy, 10 credits. Spring 2005

OPERATING INSTRUCTIONS – ENGLISH

BEFORE USING THE STIMULATOR

1. INTRODUCTION	32
2. MEDICAL BACKGROUND	33
3. PRECAUTIONARY MEASURES.....	34

USER INFORMATION

4. OVERVIEW.....	38
CONTROL BUTTONS	38
DISPLAY SYMBOLS.....	39
5. OPERATION.....	40
STEP-BY-STEP USE – PROGRAM 1 WITH HAND SWITCH.....	40
STEP-BY-STEP USE – PROGRAM 2–5.....	42
TIMER.....	43
PROGRAM PAUSE.....	43
6. REPLACEMENT OF BATTERIES.....	44
7. PROGRAM GUIDE.....	45

ADDITIONAL INFORMATION

8. ACCESSORIES	47
9. CARE INSTRUCTIONS	48
10. TROUBLESHOOTING.....	49
11. FREQUENTLY ASKED QUESTIONS (FAQ).....	50
12. TECHNICAL DATA.....	51
KEY TO THE SYMBOLS.....	52
13. ELECTRODE PLACEMENT GUIDE	SEE COVER
14. ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES.....	53
15. REFERENCES.....	57

1. INTRODUCTION

Congratulations for choosing the **CEFAR FEMINA**! Your new stimulator is especially designed for the pregnant woman of today. You can use CEFAR FEMINA before, during and after birth. CEFAR FEMINA will help you in taking care of yourself during this important period of your life.

INTENDED USE:

TENS is a well-known method based on clinical research and is widely used by medical professionals all around the world. **CEFAR FEMINA (1835)** provides a safe and natural way of pain relief during pregnancy, labour and after delivery.

More information on TENS, EMS and CEFAR products can be found on our web site: www.cefarcplex.com.

INTENDED USER:

The user of CEFAR FEMINA can be a Health Care Professional or a patient. The device should be used indoors and may be used in a healthcare facility setting or in a home environment.

2. MEDICAL BACKGROUND

TENS

TENS (Transcutaneous Electrical Nerve Stimulation) is clinically proven and used daily by physiotherapists, other caregivers and top athletes around the world. It is a safe treatment method and has, in contrast to painkillers and other pain relief methods, no side effects on you or your expected child.

Your body has its own systems for easing pain. TENS activates the pain-inhibiting mechanisms of the nervous system. Electrical impulses from electrodes, placed on the skin over or near the painful area, stimulates the nerves to block the pain signals to the brain, and the pain is not perceived. TENS also stimulates the release of endorphins, the body's natural painkillers. TENS does not always treat the cause of pain. Consult your doctor if pain persists.

EMS

Normal muscular activity is controlled by the central and peripheral nervous systems, which transmit electrical signals to the muscles. EMS (Electrical Muscle Stimulation) works similarly but uses an external source (the stimulator) with electrodes attached to the skin for transmitting electrical impulses into the body. The impulses stimulate the nerves to send signals to a specifically targeted muscle, which reacts by contracting, just as it does with normal muscular activity.

Stimulation also boosts circulation, which makes it an ideal method for tense muscles and swollen legs. The increased circulation helps the muscles to relax and improves transportation of excessive fluid from feet and legs.

3. PRECAUTIONARY MEASURES

3.1. CONTRAINDICATIONS

- Implanted electronic devices. Do not use the device on patients who have a cardiac pacemaker, implanted defibrillator, or other implanted electronic device, because this may cause electric shock, burns, electrical interference, or death.
- TENS for Undiagnosed Pain. Do not use the device as a TENS device on patients whose pain syndromes are undiagnosed.

3.2. WARNINGS

- Consult with physician. Consult with the patient's physician before using the device, because the device may cause lethal rhythm disturbances to the heart in susceptible individuals.
- Skin condition. Apply stimulation only to normal, intact, clean, healthy skin.
- Long term effects. The long-term effects of chronic electrical stimulation are unknown.
- Stimulation location. Stimulation over Neck or Mouth. Do not apply stimulation over the patient's neck (especially the carotid sinus) or the patient's mouth, because this could cause severe muscle spasms resulting in closure of the airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure.
- Stimulation across Chest. Do not apply stimulation across the patient's chest, because the introduction of electrical current into the chest may cause rhythm disturbances to the patient's heart, which could be lethal.
- Across the head. Since the effects of stimulation of the brain are unknown, stimulation should not be applied across the head, and electrodes should not be placed on opposite sides of the head.
- Stimulation over compromised skin. Do not apply stimulation over open wounds or rashes, or over swollen, red, infected, or inflamed areas or skin eruptions (e.g., phlebitis, thrombophlebitis, varicose veins).
- Stimulation near cancerous lesions. Do not apply stimulation over, or in proximity to, cancerous lesions.
- Stimulation over eyes. Do not apply stimulation directly on the eyes.
- Environment. Electronic monitoring equipment. Do not apply stimulation in the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms), which may not operate properly when the electrical stimulation device is in use.
- Bath or Shower. Do not apply stimulation when the patient is in the bath or shower. Do not apply stimulation in humid atmosphere exceeding 75% of relative humidity.
- Sleeping. Do not apply stimulation while the patient is sleeping.
- Driving or operating machinery. Do not apply stimulation while the patient is driving, operating machinery, or during any activity in which electrical stimulation or involuntary muscle contraction can put the patient at risk of injury.
- Electrosurgical equipment or defibrillators. Disconnect the stimulation electrodes before using electrosurgical equipment or defibrillators. Otherwise skin burns may be caused below the electrodes and the device might be destroyed.
- Magnetic Resonance Imaging. Do not wear electrode or the device during Magnetic Resonance Imaging (MRI) scans as this may result in metal overheating and causing skin burns in the area of the electrode.
- Flammable or explosive environment. Do not use the device in areas where there is a risk of fire or explosion, such as oxygen-rich environments, in the vicinity of flammable anaesthetics, etc.
- Power supply. Never connect stimulation cables to an external power supply as there is a risk of electric shock.
- Near other equipment. Do not use the device beside or stacked on top of any other equipment. If you must use it side by side or on top of another system, you should check that the device works properly in the chosen configuration.

- **Miscellaneous.** Electrodes for Single Patient. Do not share electrodes with other persons. All users should have individual set of electrodes to prevent undesirable skin reactions or disease transmission. Accessories. Use this device only with the leads, electrodes, and accessories recommended by the manufacturer. Use of other accessories may adversely affect the performance of the device or may result in stronger electromagnetic emissions or reduce the electromagnetic immunity of the device.
- **No Modification.** No modification of the equipment is allowed.

3.3. PRECAUTIONS

- **Supervision.** Use this device only under the continued supervision of a licensed practitioner. Electrode placement and stimulation settings should be based on the guidance of the prescribing practitioner.
- **Manufacturer.** The manufacturer does not take any responsibility for any electrode placements other than recommended.
- **Pregnancy.** The safety of electrical stimulation during pregnancy has not been established.
- **Skin irritation.** Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium (gel). The irritation may be reduced by using an alternate conductive medium or alternate electrode placement. Some patients may experience redness under the electrodes after a session. This redness usually disappears within a few hours. Advise the patient to consult the clinician if the skin redness does not disappear after a few hours. Do not start another stimulation session in the same area if the redness is still visible. Don't scratch the redness area.
- **Heart Disease.** Patients with suspected or diagnosed heart disease should follow precautions recommended by their physicians.
- **Epilepsy.** Patients with suspected or diagnosed epilepsy should follow precautions recommended by their physicians.
- **Internal Bleeding.** Use caution when the patient has a tendency to bleed internally, such as following an injury or fracture.
- **After Surgery.** Use caution following recent surgical procedures when stimulation may disrupt the patient's healing process.
- **Over uterus.** Do not place electrodes directly over the uterus or connect pairs of electrodes across the abdomen if you are pregnant. The reason is that, theoretically, the current could affect the foetus's heart (although there are no reports of it being harmful).
- **Lack of sensation.** Use caution if stimulation is applied over areas of skin that lack normal sensation. Don't apply stimulation on patient unable to express themselves.
- **Stimulation.** The stimulator should only be used with skin electrodes intended for nerve and muscle stimulation. Muscle soreness may occur after stimulation but usually disappears within a week.
- **Hot casing or batteries.** Under extreme use conditions, some parts of the casing might reach up to 109 °F (43 °C). Use caution when manipulating the batteries right after device use or when holding the device. There is no particular health risk associated with this temperature besides your comfort.
- **Children.** Keep this device out of the reach of children.
- **Electrode Size.** Do not use electrodes with an active area less than 16 cm², as there will be a risk of suffering a burn injury. Caution should always be exercised with current densities more than 2mA/cm².
- **Strangulation.** Do not wrap leadwires around your neck, and keep them out of the reach of children. Strangulation may result from entanglement in the leadwires.
- **Tripping.** Care should be used to avoid tripping on lead wires.
- **Damaged Device or Accessories.** Never use the device or any of its accessories if it is damaged (case, cables, etc.) or if the battery compartment is open as there is a risk of electric shock. Carefully inspect the lead wires and connectors prior to each use.

- Inspect Electrodes. Inspect electrodes before each use. Replace electrodes when they begin to deteriorate or lose adhesion. Poor contact between the electrodes and the patient's skin increases the risk of skin irritation or burns. Electrodes will last longer if used and stored according to instructions on electrode packaging. Attach the electrodes in such a way that their entire surface is in contact with the skin.
- Foreign bodies. Do not allow any foreign bodies (soil, water, metal, etc.) to penetrate the device and the battery compartment. Batteries. Do not carry batteries in a pocket, purse, or any other place where the terminals could become short-circuited (e.g. by way of paper clip). Intense heat could be generated and injury may result. Never open the battery cover during stimulation in order to avoid electrical shock. Remove the batteries from the device if you do not intend to use it for a prolonged period of time (more than 3 months). Extended storage of the batteries in the device might lead to batteries and device damage.
- Cable. The cable is best preserved if left attached to the stimulator between sessions. Do not jerk the cable or connection.
- Heat and Cold products. The use of heat or cold producing devices (e.g. electric heating blankets, heating pads or ice packs) may impair performance of the electrode or alter the patient's circulation/sensitivity and increase the risk of injury to the patient.
- Pulled muscles. Do not apply electrodes over pulled muscles. Using the stimulator on a previously extended muscle might further pull such muscle. The higher the stimulation intensity, the higher the risk to further overextends such muscle.
- Additional Precautions for TENS
 - TENS is not effective for pain of central origin, including headache.
 - TENS is not a substitute for pain medications and other pain management therapies.
 - TENS devices have no curative value.
 - TENS is a symptomatic treatment and, as such, suppresses the sensation of pain that would otherwise serve as a protective mechanism.
- Effectiveness of TENS is highly dependent upon patient selection by a practitioner qualified in the management of pain patients.

3.4. DANGERS

- Electrodes. Any Electrode with a minimum active area of 16 cm^2 may be used with this device. Use of an electrode with an area less than 16 cm^2 can cause burns when the unit is used at higher intensities. Consult your clinician prior to using any electrode less than 16 cm^2 .

3.5. ADVERSE REACTIONS


- Patients may experience skin irritation and burns beneath the stimulation electrodes applied to the skin.
- Patients may experience headache and other painful sensations during or following the application of electrical stimulation near the eyes and to the head and face.
- Patients should stop using the device and should consult with their physicians if they experience adverse reactions from the device.
- Precaution: Do not disconnect any stimulation cables during a session while the stimulator is switched on. Switch the stimulator off first. Always turn off the stimulator before moving or removing any electrodes during a session.
- Precaution: Do not use electrodes with an active area less than 16 cm^2 , as there will be a risk of suffering a burn injury. Caution should always be exercised with current densities more than 2 mA/cm^2 .
- Precaution: Do not apply stimulation in the vicinity of metal. Remove jewelry, body piercings, buckles or any other removable metallic product or device in the area of stimulation. Never use the electrodes contra-laterally, i.e. do not use two pins connected to the same channel on opposite segments of the body.

- Precaution: Never carry out an initial stimulation session on a person who is standing. The first five minutes of stimulation must always be performed on a person who is sitting or lying down. In rare instances, people of a nervous disposition may experience a vasovagal reaction. This reaction is connected with fear of the muscle stimulation as well as surprise at seeing one of their muscles contract without having intentionally contracted it themselves. A vasovagal reaction causes heart to slow down and blood pressure to drop, which can make you feel weak and faint. If this does occur, stop the stimulation and lie down with the legs raised until the feeling of weakness disappears (5 to 10 minutes).
- Precaution: Sudden temperature changes can cause condensation to build up inside the stimulator. To prevent this, allow it to reach room temperature before use


4. OVERVIEW/CONTROL BUTTONS

- 1. ON/OFF**
 - Turns the stimulator on and off.
 - Turns the stimulator off even when the key lock is activated. Can be used for terminating the stimulation at all times.
- 2. INCREASE (left and right channel)**
 - Increases the amplitude (intensity of stimulation). Press and hold the button to increase the amplitude continuously.
Note! Always increase the amplitude cautiously.
 - Increases the number of minutes when setting the timer (right button).
- 3. DECREASE (left and right channel)**
 - Decreases the amplitude (intensity of stimulation). Press and hold the button to decrease the amplitude continuously.
 - Deactivates the key lock.
 - Decreases the number of minutes when setting the timer (right button).
- 4. PROGRAM**
 - Selects programs 1-5.
 - Pauses an ongoing program.
- 5. TIMER**
 - Initiates the timer setting.
- 6. HAND SWITCH CONNECTION**
 - When using the hand switch, the amplitude switches between basic and intensive stimulation level in Program 1.
- 7. CABLE HOLDER**

4. OVERVIEW/DISPLAY SYMBOLS

1. BROKEN CIRCUIT

Broken circuit. The reason for a broken circuit may be too high resistance or cable breakage. See chapter TROUBLESHOOTING.

2. PAUSE

Paused program.

3. KEY LOCK

Activated key lock. The key lock is automatically activated if no key is pressed for 10 seconds. Deactivate the key lock by pressing the left or right DECREASE button.

4. BATTERY STATUS

Empty batteries. This symbol is shown when the batteries are almost empty.

5. PROGRAM NUMBER (left and right channel)

Selected program number.

6. REMAINING TIME (left and right channel)

Remaining program time in minutes and seconds. Time is flashing during timer setting.

7. AMPLITUDE BARGRAPH (left and right channel)

Selected amplitude as a bargraph.

8. AMPLITUDE LEVEL (left and right channel)

Selected amplitude in mA.

9. BASIC/INTENSIVE STIMULATION - PROGRAM 1

Indication for basic/intensive stimulation in Program 1. The upper part of the symbol is flashing during intensive stimulation and the bottom part during basic stimulation.

CONTRACTION/ACTIVE REST - PROGRAM 5

Contraction/Active Rest indication for Program 5. The upper part of the symbol is flashing during muscle contraction and the bottom part during Active Rest.


5. OPERATION

STEP-BY-STEP USE – PROGRAM 1 WITH HAND SWITCH


Program 1 is specifically designed to combat labour pain during childbirth. Two different amplitude levels are set, one for basic stimulation (between contractions) and one for intensive stimulation (during contractions). The stimulator is used together with a hand switch which lets you easily switch between the two stimulation levels.

1. INSERT THE BATTERIES

Insert the batteries (see chapter REPLACEMENT OF BATTERIES).


2. ATTACH THE ELECTRODES

A. Connect the electrodes to the cable.


B. Attach the electrodes to your body.

C. Connect the cable(s) to the CEFAR FEMINA.


3. CONNECT THE HAND SWITCH

4. SWITCH THE STIMULATOR ON

Press the **ON/OFF** button (⏻). This button can be used for terminating the stimulation at all times, even when the key lock is activated.

Always switch the stimulation off before removing electrodes from the skin.


5. SELECT PROGRAM 1

Press the **PROGRAM** button  repeatedly until P1 is shown on the display.

Note! When selecting a program the amplitude must be set to 00.0 mA for both channels.

For further information on Program 1, see chapter PROGRAM GUIDE.

6. START THE STIMULATION

A. Set the amplitude for basic stimulation (between contractions)

Press the **INCREASE** button  for each channel until you reach a comfortable level of stimulation. Press and hold the button to increase the amplitude continuously.

Note! Always increase the amplitude cautiously!


B. Set the amplitude for intensive stimulation (during contractions)

Press the hand switch button once. The basic/intensive stimulation symbol is shown on the display and its upper part is flashing. The first time you press the hand switch button, the amplitude is the same as for basic stimulation. Press the **INCREASE** button  for each channel until you reach a strong, but tolerable, stimulation.

By just pressing the hand switch button you can now quickly switch from basic stimulation between contractions to intensive stimulation when a contraction arrives. When the contraction ends, press the hand switch button again to get back to basic stimulation.

Increase the amplitude as labour progresses and pain gets more intense.

NOTE! You can preset the treatment time by using the timer function, see section TIMER.

7. STOP THE STIMULATION

Stop the stimulation by pressing the **DECREASE** button  until the amplitudes are 00.0 mA or press the **ON/OFF** button .

NOTE! If using program number 1 **without** hand switch you will only receive stimulation at one level i.e. the amplitude level set using the **INCREASE** button .

An automatically activated key lock prevents unintentional changes during treatment. The key lock is activated if no button is pressed for 10 seconds. Press any of the **DECREASE** button  to deactivate the key lock.

The last used program is saved when the stimulator is turned off and then pre-selected the next time the stimulator is turned on.


The stimulator turns off after five minutes of inactivity to spare batteries.


STEP-BY-STEP USE – PROGRAM 2-5

1. INSERT THE BATTERIES

Insert the batteries (see chapter REPLACEMENT OF BATTERIES).


2. ATTACH THE ELECTRODES

A. Connect the electrodes to the cable.

B. Attach the electrodes to your body.


C. Connect the cable(s) to the CEFAR FEMINA.

3. SWITCH THE STIMULATOR ON

Press the **ON/OFF** button . This button can be used for terminating the stimulation at all times, even when the key lock is activated.

Always switch the stimulation off before removing electrodes from the skin.


4. SELECT A PROGRAM (2-5)

Press the **PROGRAM** button  repeatedly until the program of your choice is shown on the display.

Note! When selecting programs the amplitude must be set to 00.0 mA for both channels.


For further information on the programs, see chapter PROGRAM GUIDE.

5. START THE STIMULATION

Press the **INCREASE** button  for each channel until you reach a comfortable level of stimulation (see Recommended intensity in chapter PROGRAM GUIDE). Press and hold the button to increase the amplitude continuously.

Note! Always increase the amplitude cautiously!

7. PROGRAM GUIDE

CEFAR FEMINA comprises five preset programs that are designed to alleviate pain and other discomforts that you will experience before, during and after birth. This program guide will help you to choose the right program for your varying needs during this time. It will also give you recommendations for electrode placement, body position, stimulation intensity and treatment time. To each program is at least one electrode placement number, which refers to the suggested electrode placements shown in the cover of the manual.

■ PROGRAM 1

PAIN RELIEF – CHILD BIRTH

This program is intended for pain relief during childbirth and is used together with a hand switch.

By using CEFAR FEMINA during labour and delivery you can control your pain relief yourself, you can remain mobile while using it and it does not impair your consciousness. The hand switch lets you quickly switch between basic stimulation (between contractions) and intensive stimulation (during contractions). TENS is most effective for back pain during the first stage of labour.

TENS provides safe and natural pain relief with no side effects for you or your baby. It can be used as the sole treatment or in addition to other forms of analgesics.

Recommended intensity

Basic stimulation (between contractions): Increase the amplitude gradually until you feel a distinct, yet comfortable stimulation.

Intensive stimulation (during contractions): Increase the amplitude gradually until you reach a strong but not painful stimulation.

Recommended treatment time

According to your needs.

Preset program time

No preset time is set for this program.

Electrode placement

Place the electrodes in the painful area. Try different electrode placements until the stimulation covers the whole painful area. Use 2 or 4 electrodes depending on the location of the pain and the size of the area you want to cover.

If you have pain in the front, you should place one of the electrodes in the groin area, and the other electrode (on the same cord) on the outer part of the hip or on the back. Do not place electrodes directly over the uterus or connect pairs of electrodes across the abdomen.

It is important that you change the positions of the electrodes if the painful area moves as the labour progresses.

Electrode placement number: 1, 2, 3

■ PROGRAM 2

PAIN RELIEF – PELVIS, BACK, NECK

This program can be used whenever you feel pain in pelvis, back or neck/shoulders during and after pregnancy.

These kinds of pain are common discomforts during pregnancy and are due to the increased strain put on your body and to hormonal changes that loosen your ligaments and joints.

TENS has been proven to provide effective relief from these pains and may sometimes be sufficient on its own, but can also be used as a valuable complement to other forms of treatment, for example physical exercises.

Recommended intensity

Increase the amplitude gradually until the tingling sensation is a little bit stronger than the existing pain without eliciting muscle contractions.

Recommended treatment time

Each stimulation session should last at least 30 minutes and can be performed whenever you need pain relief.

Preset program time

30 minutes.

Electrode placement

Place the electrodes on the painful area.

Electrode placement number: 3, 4, 5, 6, 7, 8, 9

■ PROGRAM 3

MASSAGE

This program can be used during pregnancy and after birth when muscles are tense.

A common reaction to stress and strain on the body is muscular tensions, often located in the back or neck/shoulders. Such tensions can increase in pregnant women and newly delivered mothers.

This massage program is inducing muscle vibrations that help tense muscles to relax and circulation in the treated area to improve.

Recommended intensity

Increase the amplitude until you get visible muscle vibrations.

Recommended treatment time

Each stimulation session should last at least 30 minutes and can be performed whenever you need massage/muscle relaxation.

Preset program time

30 minutes.

Electrode placement

Place the electrodes on the muscles that are tense.

Electrode placement number: 3, 5, 6, 7, 8, 9

■ PROGRAM 4

CIRCULATION – TIRED LEGS

This program can be used to relieve tired and swollen legs, usually occurring during the latter part of pregnancy.

These symptoms can also arise from standing or sitting in one position for a long time. The reason is often reduced circulation.

The stimulation produces a muscular pumping effect that increases the circulation in the venous and lymphatic systems. This promotes removal of waste products and excessive fluid from the treated area.

Recommended intensity

Increase the amplitude until you get visible muscle vibrations to light contractions.

Recommended treatment time

The program can be used every day.

Preset program time

30 minutes.

Electrode placement

Sit or lie down with legs elevated during stimulation.

Electrode placement number: 10

■ PROGRAM 5

GET IN SHAPE

This program can be used to get in shape after giving birth.

To restore your body to its pre-pregnancy shape, you have to do exercises that strengthen the back, abdominal and pelvic floor muscles. This program is inducing muscle contractions and will help you to regain muscle strength, vigour and elasticity. Use it as a complement to your other training.

Recommended intensity

Between the muscle contractions are periods of so called Active Rest, which means that the stimulation is active also during the rest periods. The Active Rest stimulation keeps the muscle prepared for the next contraction and reduces muscle soreness afterwards.

Contraction: Increase the amplitude gradually until you get a medium-intensity contraction.

Active Rest: Increase the amplitude gradually until you get light muscle vibrations.

Recommended treatment time

The program can be used 3-4 times a week. Also improve your fitness by walking, cycling, swimming or other strength training a couple of times per week.

Preset program time

30 minutes.

Caution! If you have had a Caesarean section you should wait 6-8 weeks before using this program, because your muscles need this time to recover.

Electrode placement number: 11, 12, 13, 14, 15, 16

9. CARE INSTRUCTIONS

Taking care of and cleaning the Cefar equipment is simple with the following instructions:

- Keep stimulator, hand switch and accessories in the original case when they are not in use. It may, however, be practical to allow the electrodes to remain on the body between treatments. In connection with stimulation, make sure that the electrodes are firmly in place.
- Self-adhesive multi-use electrodes are re-moistened if necessary with a few drops of water and kept air-tight (in a plastic bag) on protective paper when they are not in use.
- When using carbon rubber electrodes, use plenty of electrode gel and avoid drying out by applying tape around all the edges of the electrodes. Rinse the carbon rubber electrodes and the skin with water after use. Do not use detergent for the electrodes.
- Never expose the stimulator to water. Wipe it off with a damp cloth if necessary.
- Do not jerk cables or connections.
- The cables are best preserved if left attached to the stimulator between sessions.
- No other maintenance of the device is required. Service life of the device can vary depending on usage conditions. Typical service life is 7 years.
- Never service the device while in use.
- The device should be operated in temperatures between 10°C and 40°C, atmospheric pressures between 50 and 106 kPa, and relative humidity between 30% and 75%.
- The device should be transported and stored in temperatures between -40 °C and 70 °C, atmospheric pressures between 50 and 106 kPa and relative humidity between 10% and 90%.
- Contact manufacturer for assistance in setting up, using or maintaining the equipment or report events.

10. TROUBLESHOOTING

THE STIMULATION DOES NOT FEEL THE SAME AS USUAL

- Check that all settings are correct (see section STEP-BY-STEP-USE) and make sure that the electrodes are correctly placed.
- Slightly change the position of the electrodes.

THE STIMULATION FEELS UNPLEASANT

- The skin is irritated. For advice on skin care, see chapter PRECAUTIONARY MEASURE.
- The electrodes begin to lose their stickiness and do not stick properly to the skin. Moisten the adhesive surface with a few drops of water before placing on the skin.
- The electrodes are worn out and need to be replaced.
- Slightly change the position of the electrodes.

THE STIMULATION FEELS WEAK OR NOT AT ALL

- Check if the batteries need to be replaced, see chapter REPLACEMENT OF BATTERIES.
- Electrodes are too old and need replacement.

THE BROKEN CIRCUIT SYMBOL IS SHOWN ON THE DISPLAY SYMBOL

The broken circuit symbol indicates that the resistance is too high, or that a cable is broken.

- A too high resistance can be caused by a bad connection between the electrodes and your skin, or that the electrodes need to be replaced.
- A cable breakage can be checked by pressing the cable's pins against one another while increasing the amplitude for the corresponding channel to 11 mA. If the amplitude now drops to 0.0 mA and  starts flashing, the cable needs to be replaced.

Note! Never increase the amplitude above 20 mA when you check for cable breaks, since this can damage the stimulator.

THE STIMULATOR IS NOT WORKING

If the error symbol appears on the display when you start the stimulator, it means that the stimulator is broken and needs to be replaced.

Note! Do not use the stimulator – contact your Cefar dealer.

Cefar will only be responsible for service and repairs performed by Cefar or a distributor appointed by Cefar.

11. FREQUENTLY ASKED QUESTIONS (FAQ)

CAN ANYONE USE ELECTRICAL STIMULATION?

People with implanted electrical equipment for example a pacemaker and an intracardiac defibrillator must not be treated with electrical stimulation. Pregnant women should not use electrical stimulation during the first 12 weeks of the pregnancy. Read the safety precautions in this manual (PRECAUTIONARY MEASURES).

WHEN CAN I START USING TENS DURING MY PREGNANCY?

There are no medical studies showing TENS to be unsafe to use in the beginning of pregnancy. However, to be completely safe we recommend that you start using TENS after the first 12 weeks of pregnancy. Then you can use your CEFAR FEMINA regularly to get to know your stimulator and also to get in good shape for delivery.

Use the stimulator for your aching back and other pains during your pregnancy and start the stimulation at home as soon as the contractions begin. Remember that you can use TENS for afterbirth pain and post-caesarean pain as well.

WHERE DO I PLACE THE ELECTRODES?

Place the electrodes in the painful area. Try different electrode placements until the stimulation covers the whole painful area. Use 2 or 4 electrodes depending on the location of the pain and the size of the area you want to cover. A distance of 3 to 30 cm between the electrodes is recommended.

WHERE SHOULD I NOT PLACE THE ELECTRODES?

Do not place electrodes directly over the uterus or connect pairs of electrodes across the abdomen. The reason is that, theoretically, the current could affect the foetus's heart (although there are no reports of it being harmful).

If you have pain in the front, it is better to place one of the electrodes in the groin area, and the other electrode (on the same cord) on the outer part of the hip or on the back. This placement is safe, since two electrodes connected to the same cord always work together.

CAN TENS BE USED DURING CTG MONITORING?

External CTG monitoring will never interfere with TENS. However, if a scalp electrode is used, TENS may cause interference with the CTG monitor (depending on the equipment). In this case, you have to turn off the stimulator temporarily.

FOR HOW LONG DOES THE ELECTRODES LAST?

The self-adhesive electrodes last for approximately 20 to 40 occasions. The durability depends on how good the care and maintenance instructions are followed.

ARE THERE ANY SIDE EFFECTS WITH ELECTRICAL STIMULATION?

No, on the contrary to for example analgesic medicine, TENS stimulation is completely without side effects for the mother or the expected child.

ARE THERE ANY PROOFS OF TENS WORKING AS AN ALLEVIATION OF PAIN?

TENS is an established and well-documented method that has been used by the medical service for many years. Research and medical studies are carried out all the time. Read more about current studies on Cefar's home page www.cefarcompex.com

12. TECHNICAL DATA

CEFAR FEMINA is a nerve stimulator (TENS) with two non-independent channels. It has five preset programs.

Treatment with electrical stimulation requires the stimulation current to penetrate the resistance of the skin and the electrode, about 1500 ohms. CEFAR FEMINA can penetrate this resistance and maintain a current of up to 60 mA. With a change in load from 100 to 1500 ohms, the stimulation current changes less than 10% from the set value.

The stimulator operates on either two non-rechargeable 1.5 V AA batteries or two rechargeable 1.2 V AA batteries, recharged in a separate charger.

CEFAR FEMINA

Number of channels.....	2 (non-independent)
Constant current	Up to a resistance of 1500 ohm (increased load can reduce the maximum current)
Stimulation current/channel.....	0–60 mA (charge maximum: 18µC)
Waveform.....	Symmetric biphasic pulse, 100% compensated
Number of preset programs	5
Max pulse duration.....	300 µs
Max frequency.....	80 Hz
Timer.....	1 to 60 min/Off
Environment for storage, use and shipping.....	Temperature 10° C–40° C Air humidity 30%–75% Air pressure 700 hPa–1060 hPa
Power source.....	2 x 1.5 V AA non-rechargeable or 2 x 1.2 V AA rechargeable batteries
Current consumption for one channel, 80 Hz, 30 mA.....	150 mA
I r.m.s. max/channel.....	11 mA
Size	120 x 50 x 30 mm
Weight	ca. 180 g

KEY TO THE SYMBOLS

Reference number, part number

Follow instruction for uses

Type BF applied parts

Council Directive 2012/19/EU concerning Waste Electrical and Electronic Equipment (WEEE) requires not to dispose of WEEE as municipal waste. Contact your local distributor for information of the unit and accessories

Keep the device dry

Minimum and maximum temperature indications to respect

This equipment complies with all requirements of the Medical Device Directive (93/42/EEC)

Manufacturer name and address, manufacturing year

Power/Pause

14. ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
<p>CEFAR FEMINA is intended for use in the electromagnetic environment specified below. The customer or user of the CEFAR FEMINA should assure that it is used in such an environment</p>		
Emissions Tests	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Group 1	The CEFAR FEMINA uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	CEFAR FEMINA is suitable for use in any establishment, other than a private dwelling or a place connected directly to the low voltage mains supply which powers residential buildings.
Harmonic emissions IEC 61000-3-2	Not Applicable - Battery powered	
Voltage fluctuations/ emission oscillations IEC 61000-3-3	Not Applicable - Battery powered	

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

CEFAR FEMINA is intended for use in the electromagnetic environment specified below. The customer or user of the CEFAR FEMINA should assure that it is used in such an environment

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	Not Applicable - Battery powered Not Applicable - signal lines less than 3 meters	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+ 1kV differential mode (line to line) + 2kV common mode (line to ground)	Not Applicable - Battery powered	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Not Applicable - Battery powered	Mains power quality should be that of a typical commercial or hospital environment. If the user of the CEFAR FEMINA requires continued operation during power mains interruptions, it is recommended that the CEFAR FEMINA be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U_T is the AC mains voltage prior to application of the test level.			

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The CEFAR FEMINA is intended for use in the electromagnetic environment specified below. The customer or the user of the CEFAR FEMINA should assure that it is used in such an environment

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communication equipment should be used no closer to any part of the CEFAR FEMINA device, including cables, than recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> <p>$d = 1.2 \sqrt{P}$</p> <p>Radiated RF IEC 61000-4-6</p> <p>3 Vrms 80 MHz to 2.5 GHz</p> <p>3 Vrms</p> <p>$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz</p> <p>$d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CEFAR FEMINA device is used exceeds the applicable RF compliance level above, the CEFAR FEMINA device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the CEFAR FEMINA device.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 1 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the CEFAR FEMINA Electrotherapy System			
<p>The CEFAR FEMINA Electrotherapy System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the CEFAR FEMINA Electrotherapy System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the CEFAR FEMINA Electrotherapy System as recommended below, according to the maximum output power of the communications equipment.</p>			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

FCC REQUIREMENTS

Part 15 of the FCC Requirements	
<p>This device complies with Part 15 of the FCC Rules. Operation is subject to the following 2 conditions:</p>	<ul style="list-style-type: none"> • This device may not cause harmful interference • This device must accept any interference received, including the interference that may cause undesired operation.
FCC ID	T9J-RN42
Contains Transmitter Module IC	6514A-RN42

15. REFERENCES

Bundsen P, Peterson L-E, Selstam U

Pain relief in labor by transcutaneous electrical nerve stimulation. A prospective matched study.
Acta Obst Gyn Scand 1981; 60: 459-468

Grim L, Morey S

Transcutaneous electrical nerve stimulation for relief of parturition pain. A clinical report.
Physical Therapy 1985; 65: 337-340

Harrison R, Woods T, Shore M, Mathews G, Unwin A

Pain relief in labour using transcutaneous electrical nerve stimulation (TENS).
A TENS/TENS placebo controlled study in two parity groups.
Br J Obst and Gyn 1986; 93: 739-746

Kaplan B, Rabinerson D, Lurie S, Bar J, Krieser UR, Neri A

Transcutaneous electrical nerve stimulation (TENS) for adjuvant pain-relief during labor and delivery
Int Journal of Gynecology and Obstetrics 1998; 60: 251-255

Holgersson S, Johansson E

The effectiveness of TENS during labour – a survey
Research project in physiotherapy, 10 credits. Spring 2005

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15. RÉFÉRENCES

Bundsen P, Peterson L-E, Selstam U

Pain relief in labor by transcutaneous electrical nerve stimulation. A prospective matched study.
Acta Obst Gyn Scand. 1981, n° 60, pp. 459-468

Grim L, Morey S

Transcutaneous electrical nerve stimulation for relief of parturition pain. A clinical report.
Physical Therapy. 1985, n° 65, pp. 337-340

Harrison R, Woods T, Shore M, Mathews G, Unwin A

Pain relief in labour using transcutaneous electrical nerve stimulation (TENS).
A TENS/TENS placebo controlled study in two parity groups.
Br J Obst and Gyn. 1986, n° 93, pp. 739-746

Kaplan B, Rabinerson D, Lurie S, Bar J, Krieser UR, Neri A

Transcutaneous electrical nerve stimulation (TENS) for adjuvant pain-relief during labor and delivery
Int Journal of Gynecology and Obstetrics. 1998, n° 60, pp. 251-255

Holgersson S, Johansson E

The effectiveness of TENS during labour – a survey
Research project in physiotherapy, 10 crédits. Printemps 2005.

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Photo © Studio Ekvall, Helsingborg AB