CE  
2797

MD

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# 1. INTRODUCTION

## INTENDED USE:

**CEFAR REHAB X2 (2946)** is a dual channel nerve stimulator intended for both muscle rehabilitation (NMES) and pain relief (TENS). The stimulator features 27 preset programs and 3 custom programs. The channels are simultaneous, which means that a selected program applies for both channels. By using CEFAR EASY TOUCH™ the stimulation is automatically set to the correct intensity for each person.

## INTENDED USER:

The user of CEFAR REHAB X2 can be a Health Care Professional or a patient. The device should be used indoors and may be used in a healthcare facility setting or in a home environment

## 2. MEDICAL BACKGROUND

### NMES

NMES (NeuroMuscular Electrical Stimulation) is used successfully both in medical rehabilitation and as a complement to athletic training on all levels.

The goal of electrical muscle stimulation is to achieve contractions or vibrations in the muscles. Normal muscular activity is controlled by the central and peripheral nervous systems, which transmit electrical signals to the muscles. NMES works similarly but uses an external source (the stimulator) with electrodes attached to the skin for transmitting electrical impulses into the body. The impulses stimulate the nerves to send signals to a specifically targeted muscle, which reacts by contracting, just as it does with normal muscular activity.

Electrical muscle stimulation is suitable for all the muscles in the body. It can be used to strengthen muscles weakened by surgery, a fracture, etc., and improve mobility. It is also an excellent tool for stroke rehabilitation, helping patients in handgrip and gait training.

Electrical muscle stimulation for rehabilitation purposes should be tried out individually by a physiotherapist or other caregiver for the best results.

### TENS

TENS (Transcutaneous Electrical Nerve Stimulation) gives good results in acute and chronic pain conditions of many kinds. It is clinically proven and used daily by physiotherapists, other caregivers and top athletes around the world.

High-frequency TENS activates the pain-inhibiting mechanisms of the nervous system. Electrical impulses from electrodes, placed on the skin over or near the painful area, stimulate the nerves to block the pain signals to the brain, and the pain is not perceived. Low-frequency TENS stimulates the release of endorphins, the body's natural painkillers.

TENS is a safe treatment method and has, in contrast to drugs and other pain relief methods, no side effects. It may be sufficient as the only treatment form, but it is also a valuable complement to other pharmacological and/or physical treatments. TENS does not always treat the cause of pain. Consult your doctor if pain persists.

## 3. PRECAUTIONARY MEASURES

### 3.1. CONTRAINDICATIONS

- Implanted electronic devices. Do not use the device on patients who have a cardiac pacemaker, implanted defibrillator, or other implanted electronic device, because this may cause electric shock, burns, electrical interference, or death.
- TENS for Undiagnosed Pain. Do not use the device as a TENS device on patients whose pain syndromes are undiagnosed.

### 3.2. WARNINGS

- Consult with physician. Consult with the patient's physician before using the device, because the device may cause lethal rhythm disturbances to the heart in susceptible individuals.
- Skin condition. Apply stimulation only to normal, intact, clean, healthy skin.
- Long term effects. The long-term effects of chronic electrical stimulation are unknown.
- Stimulation location. Stimulation over Neck or Mouth. Do not apply stimulation over the patient's neck (especially the carotid sinus) or the patient's mouth, because this could cause severe muscle spasms resulting in closure of the airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure.
- Stimulation across Chest. Do not apply stimulation across the patient's chest, because the introduction of electrical current into the chest may cause rhythm disturbances to the patient's heart, which could be lethal.
- Across the head. Since the effects of stimulation of the brain are unknown, stimulation should not be applied across the head, and electrodes should not be placed on opposite sides of the head.
- Stimulation over compromised skin. Do not apply stimulation over open wounds or rashes, or over swollen, red, infected, or inflamed areas or skin eruptions (e.g., phlebitis, thrombophlebitis, varicose veins).
- Stimulation near cancerous lesions. Do not apply stimulation over, or in proximity to, cancerous lesions.
- Stimulation over eyes. Do not apply stimulation directly on the eyes.
- Environment. Electronic monitoring equipment. Do not apply stimulation in the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms), which may not operate properly when the electrical stimulation device is in use.
- Bath or Shower. Do not apply stimulation when the patient is in the bath or shower. Do not apply stimulation in humid atmosphere exceeding 75% of relative humidity.
- Sleeping. Do not apply stimulation while the patient is sleeping.
- Driving or operating machinery. Do not apply stimulation while the patient is driving, operating machinery, or during any activity in which electrical stimulation or involuntary muscle contraction can put the patient at risk of injury.
- Electrosurgical equipment or defibrillators. Disconnect the stimulation electrodes before using electrosurgical equipment or defibrillators. Otherwise skin burns may be caused below the electrodes and the device might be destroyed.
- Magnetic Resonance Imaging. Do not wear electrode or the device during Magnetic Resonance Imaging (MRI) scans as this may result in metal overheating and causing skin burns in the area of the electrode.
- Flammable or explosive environment. Do not use the device in areas where there is a risk of fire or explosion, such as oxygen-rich environments, in the vicinity of flammable anaesthetics, etc.
- Power supply. Never connect stimulation cables to an external power supply as there is a risk of electric shock.
- Near other equipment. Do not use the device beside or stacked on top of any other equipment. If you must use it side by side or on top of another system, you should check that the device works properly in the chosen configuration.

- **Miscellaneous. Electrodes for Single Patient.** Do not share electrodes with other persons. All users should have individual set of electrodes to prevent undesirable skin reactions or disease transmission.
- **Accessories.** Use this device only with the leads, electrodes, and accessories recommended by the manufacturer. Use of other accessories may adversely affect the performance of the device or may result in stronger electromagnetic emissions or reduce the electromagnetic immunity of the device.
- **No Modification.** No modification of the equipment is allowed.

### 3.3. PRECAUTIONS

- **Supervision.** Use this device only under the continued supervision of a licensed practitioner. Electrode placement and stimulation settings should be based on the guidance of the prescribing practitioner.
- **Manufacturer.** The manufacturer does not take any responsibility for any electrode placements other than recommended.
- **Pregnancy.** The safety of electrical stimulation during pregnancy has not been established.
- **Skin irritation.** Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium (gel). The irritation may be reduced by using an alternate conductive medium or alternate electrode placement. Some patients may experience redness under the electrodes after a session. This redness usually disappears within a few hours. Advise the patient to consult the clinician if the skin redness does not disappear after a few hours. Do not start another stimulation session in the same area if the redness is still visible. Don't scratch the redness area.
- **Heart Disease.** Patients with suspected or diagnosed heart disease should follow precautions recommended by their physicians.
- **Epilepsy.** Patients with suspected or diagnosed epilepsy should follow precautions recommended by their physicians.
- **Internal Bleeding.** Use caution when the patient has a tendency to bleed internally, such as following an injury or fracture.
- **After Surgery.** Use caution following recent surgical procedures when stimulation may disrupt the patient's healing process.
- **Over uterus.** Do not place electrodes directly over the uterus or connect pairs of electrodes across the abdomen if you are pregnant. The reason is that, theoretically, the current could affect the foetus's heart (although there are no reports of it being harmful).
- **Lack of sensation.** Use caution if stimulation is applied over areas of skin that lack normal sensation. Don't apply stimulation on patient unable to express themselves.
- **Stimulation.** The stimulator should only be used with skin electrodes intended for nerve and muscle stimulation. Muscle soreness may occur after stimulation but usually disappears within a week.
- **Hot casing or batteries.** Under extreme use conditions, some parts of the casing might reach up to 109 °F (43 °C). Use caution when manipulating the batteries right after device use or when holding the device. There is no particular health risk associated with this temperature besides your comfort.
- **Children.** Keep this device out of the reach of children.
- **Electrode Size.** Do not use electrodes with an active area less than 16 cm<sup>2</sup>, as there will be a risk of suffering a burn injury. Caution should always be exercised with current densities more than 2mA/ cm<sup>2</sup>.
- **Strangulation.** Do not wrap leadwires around your neck, and keep them out of the reach of children. Strangulation may result from entanglement in the leadwires.
- **Tripping.** Care should be used to avoid tripping on lead wires.
- **Damaged Device or Accessories.** Never use the device or any of its accessories if it is damaged (case, cables, etc.) or if the battery compartment is open as there is a risk of electric shock. Carefully inspect the lead wires and connectors prior to each use.

- **Inspect Electrodes.** Inspect electrodes before each use. Replace electrodes when they begin to deteriorate or lose adhesion. Poor contact between the electrodes and the patient's skin increases the risk of skin irritation or burns. Electrodes will last longer if used and stored according to instructions on electrode packaging. Attach the electrodes in such a way that their entire surface is in contact with the skin.
- **Foreign bodies.** Do not allow any foreign bodies (soil, water, metal, etc.) to penetrate the device and the battery compartment.
- **Batteries.** Do not carry batteries in a pocket, purse, or any other place where the terminals could become short-circuited (e.g. by way of paper clip). Intense heat could be generated and injury may result. Never open the battery cover during stimulation in order to avoid electrical shock. Remove the batteries from the device if you do not intend to use it for a prolonged period of time (more than 3 months). Extended storage of the batteries in the device might lead to batteries and device damage.
- **Cable.** The cable is best preserved if left attached to the stimulator between sessions. Do not jerk the cable or connection.
- **Heat and Cold products.** The use of heat or cold producing devices (e.g. electric heating blankets, heating pads or ice packs) may impair performance of the electrode or alter the patient's circulation/sensitivity and increase the risk of injury to the patient.
- **Pulled muscles.** Do not apply electrodes over pulled muscles. Using the stimulator on a previously extended muscle might further pull such muscle. The higher the stimulation intensity, the higher the risk to further overextends such muscle.
- **Additional Precautions for TENS**
  - TENS is not effective for pain of central origin, including headache.
  - TENS is not a substitute for pain medications and other pain management therapies.
  - TENS devices have no curative value.
  - TENS is a symptomatic treatment and, as such, suppresses the sensation of pain that would otherwise serve as a protective mechanism.
- Effectiveness of TENS is highly dependent upon patient selection by a practitioner qualified in the management of pain patients.

### 3.4. DANGERS

- **Electrodes.** Any Electrode with a minimum active area of 16 cm<sup>2</sup> may be used with this device. Use of an electrode with an area less than 16 cm<sup>2</sup> can cause burns when the unit is used at higher intensities. Consult your clinician prior to using any electrode less than 16 cm<sup>2</sup>.

### 3.5. ADVERSE REACTIONS

- Patients may experience skin irritation and burns beneath the stimulation electrodes applied to the skin.
- Patients may experience headache and other painful sensations during or following the application of electrical stimulation near the eyes and to the head and face.
- Patients should stop using the device and should consult with their physicians if they experience adverse reactions from the device.
- **Precaution:** Do not disconnect any stimulation cables during a session while the stimulator is switched on. Switch the stimulator off first. Always turn off the stimulator before moving or removing any electrodes during a session.
- **Precaution:** Do not use electrodes with an active area less than 16 cm<sup>2</sup>, as there will be a risk of suffering a burn injury. Caution should always be exercised with current densities more than 2mA/cm<sup>2</sup>.
- **Precaution:** Do not apply stimulation in the vicinity of metal. Remove jewelry, body piercings, buckles or any other removable metallic product or device in the area of stimulation. Never use the electrodes contra-laterally, i.e. do not use two pins connected to the same channel on opposite segments of the body.

- Precaution: Never carry out an initial stimulation session on a person who is standing. The first five minutes of stimulation must always be performed on a person who is sitting or lying down. In rare instances, people of a nervous disposition may experience a vasovagal reaction. This reaction is connected with fear of the muscle stimulation as well as surprise at seeing one of their muscles contract without having intentionally contracted it themselves. A vasovagal reaction causes heart to slow down and blood pressure to drop, which can make you feel weak and faint. If this does occur, stop the stimulation and lie down with the legs raised until the feeling of weakness disappears (5 to 10 minutes).
- Precaution: Sudden temperature changes can cause condensation to build up inside the stimulator. To prevent this, allow it to reach room temperature before use.

## 4. OVERVIEW

### CONTROL BUTTONS

#### 1. ON/OFF

- Turns the stimulator on and off.
- Turns the stimulator off even when the key lock is activated.
- Can be used for terminating the stimulation at all times.

#### 2. INCREASE (left and right channel)

- Increases the amplitude (intensity of stimulation). Press and hold the button to increase the amplitude continuously.
- **Note!** Always increase the amplitude cautiously.
- Terminates AUTO stimulation mode.
- Increases the number of minutes when setting the timer (right button).
- Used when turning the program lock on and off (left button)\*.
- Scrolls through choices in programming mode\*.

#### 3. DECREASE (left and right channel)


- Decreases the amplitude (intensity of stimulation). Press and hold the button to decrease the amplitude continuously.
- Deactivates the key lock.
- Terminates AUTO stimulation mode.
- Decreases the number of minutes when setting the timer (right button).
- Used when turning the program lock on and off (left button)\*.
- Scrolls through choices in programming mode\*.

#### 4. AUTO (left and right channel)

- Starts the AUTO test mode.
- Confirms the selected amplitude level in AUTO test mode.

#### 5. PROGRAM

- Selects a program (P1-P30). Use P+ to step forward through the programs and P- to step backward.
- Pauses an ongoing program.
- Used when activating/deactivating the program lock (P+)\*.

#### 6. TIMER

- Initiates the timer setting.

#### 7. PROGRAMMING/CONFIRMATION (S=SET)

- Turns the stimulator into programming mode for custom programs P28-P30 when pressed for 2 seconds.
- Confirms settings in programming mode.

#### 8. HAND SWITCH CONNECTION

- By using the hand switch, the contractions in programs with intermittent stimulation can be manually controlled.

#### 9. CABLE HOLDER

\*Professional use


# 5. OPERATION

## STEP-BY-STEP USE


CEFAR REHAB X2 can be used for both TENS and NMES treatment. Use only those treatment programs that your health care provider has determined are the most effective for your needs. The stimulator has two simultaneous channels, which means that both channels stimulate with the same program. A hand switch is provided to make it easy to manually control muscle contractions in NMES programs.

### 1. INSERT THE BATTERIES


Insert the batteries (see chapter *REPLACEMENT OF BATTERIES*).


### 2. ATTACH THE ELECTRODES

A. Connect the electrodes to the cable.


B. Attach the electrodes to your body.

C. Connect the cable to the CEFAR REHAB X2.

### 3. CONNECT THE HAND SWITCH


**Note!** This is only applicable for programs P9–P16 and P19–P23.



### 4. SWITCH THE STIMULATOR ON

Press the **ON/OFF** button . This button can be used for terminating the stimulation at all times, even when the key lock is activated.

Always switch the stimulation off before removing electrodes from the skin.


### 5. SELECT A PROGRAM (P1-P30)

Press the **PROGRAM** button  or  to step forward or backward until the program of your choice is shown on the display.

**Note!** When selecting a program the amplitude must be set to 00.0 mA for both channels.

For further information on the programs, see chapter *PROGRAMS*.

### 6. START THE STIMULATION

#### TENS


Press the **INCREASE** button  for each channel until you reach a comfortable level of stimulation. Press and hold the button to increase the amplitude continuously.

#### NMES

Press the **INCREASE** button  for each channel until you get visible muscle contractions. Press and hold the button to increase the amplitude continuously.

**Note!** Both channels must be used for programs P5 and P17-P18.

**Note!** Always increase the amplitude cautiously!

See section SPECIAL INSTRUCTIONS for the following cases:

- Mixed frequency stimulation
- Intermittent stimulation + hand switch
- Intermittent stimulation with Active Rest

An automatically activated key lock prevents unintentional changes during treatment. The key lock is activated if no button is pressed for 20 seconds. Press any of the **DECREASE** buttons  to deactivate the key lock.

The treatment time is preset, but can be changed by using the timer function, see section *TIMER*.

### 7. STOP THE STIMULATION

You can stop the stimulation before the program time has elapsed by pressing the **DECREASE** button  until the amplitudes are 00.0 mA or by pressing the **ON/OFF** button .

When the program is finished, remaining time "00:00" is flashing on the display.

The last used program is saved when the stimulator is turned off and then pre-selected the next time the stimulator is turned on.

The stimulator turns off after 5 minutes of inactivity to spare batteries.

## SPECIAL INSTRUCTIONS

### Mixed frequency stimulation (P4)

Mixed frequency stimulation is a combination of high (80 Hz) and low (2 Hz) frequency stimulation with separate amplitude settings.


1. First set the amplitude for 80 Hz stimulation when the upper part of the high/low frequency symbol  is flashing
2. When the bottom part of the symbol flashes for the first time, the amplitude for 2 Hz stimulation will start at the same level as is set for 80 Hz. Adjust the amplitude until you get visible muscle contractions.

The period between the frequency changes is preset to 3 seconds.

See chapter *FREQUENTLY ASKED QUESTIONS* for further information on mixed frequency stimulation.

### Intermittent stimulation + hand switch (P9–P23)

Programs 9–23 are intermittent stimulation programs with rest time between the muscle contractions (work time), see figure. The contraction gradually increases during ramp-up time and reaches maximum contraction during work time. During ramp-down time, the stimulation gradually decreases until the rest time starts. The rest time can be without or with stimulation (Active Rest). In intermittent stimulation programs the work/rest symbol  is shown on the display. The upper part of the symbol is flashing during work time and the bottom part during rest time.

Programs P9–P16 and P19–P23 can be used together with the hand switch to manually control the duration of the work time and rest time. Press the hand switch button during rest time for ramp-up time and contraction. The duration of the work time is preset but can be shortened by pressing the hand switch button. The rest time lasts until the hand switch button is pressed again.

See chapter *FREQUENTLY ASKED QUESTIONS* for further information on using the hand switch.

### Intermittent stimulation with Active Rest (P19–P20)

In programs with Active Rest, the stimulation is active also during the rest time, which means that you have to perform two amplitude settings for each channel.

**Amplitude setting for contractions (work time):** When the upper part of the work/rest symbol  is flashing, increase the amplitude gradually until you get muscle contractions without pain.

**Amplitude setting for Active Rest:** When the bottom part of the work/rest symbol  is flashing, increase the amplitude gradually until you get muscle vibrations.

See chapter *FREQUENTLY ASKED QUESTIONS* for further information on Active Rest.


# 6. REPLACEMENT OF BATTERIES

A battery symbol is shown on the display when the batteries are almost empty . As long as the stimulator is working normally you can continue the treatment. When stimulation feels weaker than usual or the stimulator turns off, it is time to replace the batteries.

If the stimulator is not used for some time (approximately 3 months), the batteries should be removed from the stimulator.

The stimulator operates on either two non-rechargeable 1.5 V AA batteries or two rechargeable 1.2 V AA batteries, recharged in a separate charger.

## REPLACING THE BATTERIES

1. Turn off the stimulator.
2. Locate the battery compartment on the back of the stimulator.
3. Remove the battery cover by pressing the mark at the top firmly and sliding the cover downwards.
4. Remove the batteries.
5. Insert the new batteries correctly according to the polarity markings (+ and -) inside the battery compartment, see picture.
6. Replace the battery cover.
7. Dispose of the exhausted batteries in accordance with local and national regulations.

**NOTE!** Non-rechargeable batteries may not be charged due to the risk of explosion.

## 7. PROGRAMS

### PRESET PROGRAMS – TENS

CEFAR REHAB X2 has six preset TENS programs for pain relief. The program overview below shows the stimulation mode and parameter settings for each program. For information on treating various pain conditions with TENS, see the clinical guide last in this manual. For further information on programs P4–P5, see chapter *FREQUENTLY ASKED QUESTIONS*.

#### **P1 – CONVENTIONAL TENS**

High-frequency stimulation that causes a tingling sensation. This is often the first choice for both acute and long-term pain.

Frequency	80 Hz
Pulse duration	180 $\mu$ s
Total program time	-- min.

#### **P2 – BURST TENS**

Low-frequency stimulation that produces visible muscle contractions. The pain relief sets in slower but lasts longer than conventional TENS.

Frequency	2 Hz
Pulse duration	180 $\mu$ s
Total program time	-- min.

#### **P3 – MODULATED PULSE DURATION TENS**

A type of high-frequency stimulation where the pulse duration varies continuously. This can cause an undulating sensation, which may be more pleasant than constant pulse duration.

Frequency	80 Hz
Pulse duration 1	70 $\mu$ s
Pulse duration 2	180 $\mu$ s
Modulation time	2 sec.
Total program time	-- min.

#### **P4 – MIXED FREQUENCY TENS**

A combination of high- and low-frequency stimulation, which can lead to a more effective pain relief treatment.

Frequency 1	80 Hz
Frequency 2 (burst)	2 Hz
Pulse duration	180 $\mu$ s
Modulation time	3/3 sec.
Total program time	-- min.

#### **P5 – FLOW TENS**

An alternating stimulation causing a massaging and pumping effect that can be used for pain relief and increased circulation in the stimulated area. Always use 2 channels and 4 electrodes.

Frequency	80 Hz
Pulse duration 1	70 $\mu$ s
Pulse duration 2	180 $\mu$ s
Modulation time	2 sec.
Total program time	-- min.

#### **P6 – GUIDE TO INTRAMUSCULAR INJECTIONS**

A program used to increase the precision of intramuscular botox (botulinum toxin) injections. These injections are used for spasticity, dystonia, facial hemispasm and bruxism.

Frequency	1 Hz
Pulse duration	50 $\mu$ s
Total program time	-- min.

**PRESET PROGRAMS – NMES**

CEFAR REHAB X2 has 21 preset NMES programs, covering warm up, rehabilitation, strengthening, recovery and massage. In the clinical guide, last in this manual, is a list of indications with suggested treatment programs and electrode placements for each indication. P9–P16 and P19–P23 can be used together with the hand switch, which makes it easy to manually control muscle contractions.

**P7 – WARM UP, SHORT**

Frequency 1	2 Hz
Frequency 2	8 Hz
Pulse duration	300 µs
Modulation time	10 sec.
Total program time	-- min.

**P8 – WARM UP, LONG**

Frequency 1	2 Hz
Frequency 2	8 Hz
Pulse duration	300 µs
Modulation time	10 sec.
Total program time	-- min.

**P9 – INTERMITTENT STIMULATION**  
**Smaller muscle groups**

Frequency	50 Hz
Pulse duration	200 µs
Work time	4 sec.
Rest time	4 sec.
Ramp-up time	2 sec.
Ramp-down time	2 sec.
Total program time	-- min.

**P10 – INTERMITTENT STIMULATION**  
**Larger muscle groups**

Frequency	50 Hz
Pulse duration	400 µs
Work time	4 sec.
Rest time	4 sec.
Ramp-up time	2 sec.
Ramp-down time	2 sec.
Total program time	-- min.

**P11 – INTERMITTENT STIMULATION**  
**Smaller muscle groups**

Frequency	65 Hz
Pulse duration	200 µs
Work time	6 sec.
Rest time	8 sec.
Ramp-up time	2 sec.
Ramp-down time	2 sec.
Total program time	-- min.

**P12 – INTERMITTENT STIMULATION**  
**Larger muscle groups**

Frequency	65 Hz
Pulse duration	400 µs
Work time	6 sec.
Rest time	8 sec.
Ramp-up time	2 sec.
Ramp-down time	2 sec.
Total program time	-- min.

**P13 – INTERMITTENT STIMULATION**  
**Smaller muscle groups**

Frequency	50 Hz
Pulse duration	200 µs
Work time	10 sec.
Rest time	10 sec.
Ramp-up time	1 sec.
Ramp-down time	2 sec.
Total program time	-- min.

**P14 – INTERMITTENT STIMULATION****Larger muscle groups**

Frequency	50 Hz
Pulse duration	400 µs
Work time	10 sec.
Rest time	10 sec.
Ramp-up time	1 sec.
Ramp-down time	2 sec.
Total program time	-- min.

**P15 – INTERMITTENT STIMULATION****Neurology, smaller muscle groups**

Frequency	40 Hz
Pulse duration	200 µs
Work time	5 sec.
Rest time	15 sec.
Ramp-up time	4 sec.
Ramp-down time	2 sec.
Total program time	-- min.

**P16 – INTERMITTENT STIMULATION****Neurology, larger muscle groups**

Frequency	40 Hz
Pulse duration	400 µs
Work time	5 sec.
Rest time	15 sec.
Ramp-up time	4 sec.
Ramp-down time	2 sec.
Total program time	-- min.

**P17 – ALTERNATING STIMULATION****Smaller muscle groups**

Frequency	50 Hz
Pulse duration	200 µs
Work time	4 sec.
Rest time	6 sec.
Ramp-up time	1 sec.
Ramp-down time	1 sec.
Total program time	-- min.

**P18 – ALTERNATING STIMULATION****Larger muscle groups**

Frequency	50 Hz
Pulse duration	400 µs
Work time	4 sec.
Rest time	6 sec.
Ramp-up time	1 sec.
Ramp-down time	1 sec.
Total program time	-- min.

**P19 – STIMULATION WITH ACTIVE REST****Smaller muscle groups**

Frequency – work	50 Hz
Frequency – rest	8 Hz
Pulse duration	200 µs
Work time	10 sec.
Rest time	10 sec.
Ramp-up time – work	2 sec.
Ramp-down time – work	2 sec.
Ramp-up time – rest	1 sec.
Ramp-down time – rest	1 sec.
Total program time	-- min.

**P20 – STIMULATION WITH ACTIVE REST****Larger muscle groups**

Frequency – work	50 Hz
Frequency – rest	8 Hz
Pulse duration	400 µs
Work time	10 sec.
Rest time	10 sec.
Ramp-up time – work	2 sec.
Ramp-down time – work	2 sec.
Ramp-up time – rest	1 sec.
Ramp-down time – rest	1 sec.
Total program time	-- min.

**P21 – STRENGTHENING**

**Upper extremities**

Frequency 1	25 Hz
Frequency 2	45 Hz
Pulse duration	200 µs
Work time	5 sec.
Rest time	8 sec.
Ramp-up time	2 sec.
Ramp-down time	1 sec.
Total program time	-- min.

**P22 – STRENGTHENING**

**Lower extremities**

Frequency 1	45 Hz
Frequency 2	65 Hz
Pulse duration	400 µs
Work time	5 sec.
Rest time	8 sec.
Ramp-up time	2 sec.
Ramp-down time	1 sec.
Total program time	-- min.

**P23 – STRENGTHENING**

**Back and trunk**

Frequency 1	30 Hz
Frequency 2	50 Hz
Pulse duration	300 µs
Work time	5 sec.
Rest time	8 sec.
Ramp-up time	2 sec.
Ramp-down time	1 sec.
Total program time	-- min.

**P24 – RECOVERY**

**Smaller muscle groups**

Frequency	4 Hz
Pulse duration	200 µs
Total program time	-- min.

**P25 – RECOVERY**

**Larger muscle groups**

Frequency	4 Hz
Pulse duration	400 µs
Total program time	-- min.

**P26 – MASSAGE**

**Smaller muscle groups**

Frequency 1	5 Hz
Frequency 2	15 Hz
Pulse duration	200 µs
Modulation time	5 sec.
Total program time	-- min.

**P27 – MASSAGE**

**Larger muscle groups**



Frequency 1	5 Hz
Frequency 2	15 Hz
Pulse duration	400 µs
Modulation time	5 sec.
Total program time	-- min.

## CUSTOM PROGRAMS

With the CEFAR REHAB X2 it is possible to create and store three custom programs (P28–P30) for patient-specific treatment. To create a custom program, follow the programming procedure below. To use a custom program, follow the instructions in the section STEP-BY-STEP USE.

### Programming

1. Press the **ON/OFF** button,  to turn the stimulator on.
2. Press the **PROGRAM** button  or  to step forward or backward until Program 28, 29 or 30 is shown on the display. Select one of these programs.
3. Press the **PROGRAMMING/CONFIRMATION** button  for 2 seconds to enter the programming mode. (See the programming chart for the following steps.)

4. **Level 1** (shown in the upper right corner of the display):

The first step in the programming procedure is to choose between:

- Continuous (**C**) stimulation
- Intermittent (**I**) stimulation

Press the **INCREASE** (or **DECREASE**) button   to toggle between **C** and **I**. Confirm your choice by pressing the **PROGRAMMING/CONFIRMATION** button .

You are moved to the next stage in the programming procedure.

5. **Level 2**

For Continuous stimulation there are no choices at Level 2. You are moved directly to Level 3.

If you chose Intermittent stimulation on Level 1, you now have two choices:

- Simultaneous (**SI**) stimulation
- Alternating (**Alt**) stimulation

Press the **INCREASE** (or **DECREASE**) button   to toggle between **SI** and **Alt**. Confirm your choice by pressing the **PROGRAMMING/CONFIRMATION** button .

You are moved to the next stage in the programming procedure.

6. **Level 3**

The available stimulation types on this level depend on your previous choices.

Press the **INCREASE** (or **DECREASE**) button   to see the different stimulation types alternate on the display. The programming chart shows the parameter setting possibilities for each stimulation type.

Continuous stimulation types:

- Conventional (**C**)
- Burst (**B**)
- Pulse width/duration modulation (**PWM**)
- Frequency modulation (**FM**)


Intermittent/Simultaneous stimulation types:


- Conventional (**C**)
- Frequency modulation (**FM**)

For Intermittent/Alternating stimulation your only choice is Conventional stimulation, which is not

shown on the display. You are moved directly to the next stage in the programming procedure (see step 7).

Press the **PROGRAMMING/CONFIRMATION** button  to confirm your choice of stimulation type and you are moved to the next stage in the programming procedure. In the following steps you will set the parameters for the selected stimulation type.


7. A preset parameter value is shown on the screen, but you can change this value by pressing the **INCREASE** (or **DECREASE**) button  . The valid range for the parameter is shown in the programming chart.
8. Press the **PROGRAMMING/CONFIRMATION** button  to confirm your setting. You are now moved to the next stage (if any) in the programming procedure.
9. Repeat step 7 and 8 until the parameter value in the last stage is set.
10. The programming procedure is finished and you will automatically exit the programming mode. The new program is now selected and ready for use. The program is also stored in the stimulator for future use.

A custom program can be changed by redoing the programming.


## Programming chart

24

\* Frequency between 1 and 5 Hz ⇔ Burst stimulation  
Other frequency ⇔ Conventional stimulation

\*\*\* Rest time must be greater than  
Work time + Ramp-up time work + Ramp-down time work

\*\*\*\* 0.5 s is shown as 05 on the display.

## 8. ACCESSORIES



CEFAR REHAB X2 comes with a neck strap and a belt clip, allowing you to wear the stimulator around your neck or on your belt and have your hands free during treatment.

The electrodes will eventually wear out and have to be replaced. It is recommended to replace the electrodes after approximately 20–40 times of usage.

The hand switch included can be used in some NMES programs for full manual control of the time of contractions.

The cables are best preserved if left attached to the stimulator between sessions.

For purchase information, contact your Cefar dealer or visit [www.cefar.se](http://www.cefar.se).

## 9. CARE INSTRUCTIONS

Taking care of and cleaning the Cefar equipment is simple with the following instructions:

- Keep stimulator and accessories in the original case when they are not in use. It may, however, be practical to allow the electrodes to remain on the body between treatments. Carbon rubber electrodes can generally remain for 2–3 hours without the electrode gel drying out (does not apply to adhesive gel). They must then be taken off, washed, and dried before being applied again. This is especially important for persons with sensitive skin. In connection with stimulation, make sure that the electrodes are firmly in place.
- When using carbon rubber electrodes, use plenty of electrode gel and avoid drying out by applying tape around all the edges of the electrodes. Rinse the carbon rubber electrodes and the skin with water after use. Do not use detergent for the electrodes.
- Self-adhesive multi-use electrodes are re-moistened if necessary with a few drops of water and kept air-tight (in a plastic bag) on protective paper when they are not in use.
- Never expose the stimulator to water. Wipe it off with a damp cloth if necessary.
- Do not jerk cables or connections.
- The cables are best preserved if left attached to the stimulator between sessions.
- No other maintenance of the device is required. Service life of the device can vary depending on Usage conditions. Typical service life is 7 years.
- Never service the device while in use.
- The device should be operated in temperatures between 10°C and 40°C, atmospheric pressures between 50 and 106 kPa, and relative humidity between 30% and 75%.
- The device should be transported and stored in temperatures between -40 °C and 70 °C, atmospheric pressures between 50 and 106 kPa and relative humidity between 10% and 90%.
- Contact manufacturer for assistance in setting up, using or maintaining the equipment or report events.

## 11. FREQUENTLY ASKED QUESTIONS (FAQ)

### CAN ANYONE USE ELECTRICAL STIMULATION?

People with implanted electrical equipment, for example a pacemaker or an intracardiac defibrillator, must not be treated with electrical stimulation. Pregnant women should not use electrical stimulation during the first 12 weeks of the pregnancy. Read the safety precautions in this manual (PRECAUTIONARY MEASURES).

### WHEN SHOULD I USE MIXED-FREQUENCY STIMULATION?

With mixed-frequency stimulation both muscle nerves (2 Hz) and sensory nerves (80 Hz) are stimulated. The stimulator switches between these two frequencies every three seconds, providing the benefits of both high-frequency stimulation (fast pain relief) and low-frequency stimulation (longer lasting but slower pain relief). This stimulation type can lead to more effective treatment of long-term pain conditions.

### WHAT KIND OF PROGRAM IS P5 (FLOW TENS)?

Cefar Flow TENS is a new comfortable way of stimulation requiring four electrodes. Use it on large areas to obtain pain relief together with a massage/muscle relaxation. The alternating stimulation gives a pumping effect that increases the circulation in the area.

### WHAT IS THE ADVANTAGE OF USING THE HAND SWITCH?

It helps you to get an effective and individual stimulation by allowing you to manually control the time of contractions in intermittent stimulation programs. By pressing the hand switch button during contraction, the stimulation gradually decreases until the rest time starts. If you do not stop the stimulation manually during contraction, it will continue the whole preset work time. With the hand switch connected, the rest time will last until you press the hand switch button. Without hand switch the rest time is preset and depending on the program you use.

**Note!** The preset work time can be shortened but not extended. It is therefore preferable to use a program with a longer work time when using the hand switch.

**Note!** When using CEFAR EASY TOUCH™ the AUTO stimulation mode must have started before you press the hand switch button.

**Note!** The hand switch does not work during alternating stimulation (P17-18).

### WHAT DOES ACTIVE REST MEAN?

It means that low frequency stimulation is active during rest time, causing muscle vibrations to maintain circulation. The Active Rest stimulation helps eliminate lactic acid and waste products, thereby reducing muscle soreness afterwards and keeping the muscle prepared for the next contraction.

**Note!** The amplitude level must be set for both contractions and Active Rest.

### HOW LONG WILL THE ELECTRODES LAST?

The self-adhesive electrodes last for approximately 20 to 40 occasions. The durability depends on how good the care and maintenance instructions are followed.

**WHICH DISTANCE SHOULD I HAVE BETWEEN THE ELECTRODES?**

It is recommended to have a distance of 3 to 30 cm between the electrodes.

**HOW DO I FIND THE OPTIMAL POSITION OF THE ELECTRODES FOR NMES?**

Use carbon rubber electrodes and gel. Slide the electrodes slowly over the muscle while stimulating at 2 Hz. The optimal position for the electrodes is where the strongest motor response occur.

**FOR HOW LONG CAN I STIMULATE?**

**TENS (80 Hz):** Can be used without an upper limit, but at least 30 min at each occasion.

**TENS (2 Hz):** Can cause sore muscles but normally 20–45 min three times a day is recommended.

**NMES:** Depending on the muscle's status and where the patient is in the rehabilitation process, treatment can last from 5 to 60 minutes and be repeated from three times a week to twice a day. Remember that the patient may develop sore muscles after NMES treatment.

## 12. TECHNICAL DATA

**CEFAR REHAB X2** is a dual channel stimulator intended for both muscle rehabilitation (NMES) and pain relief (TENS). The stimulator features 27 preset programs and 3 custom programs.

Treatment with electrical stimulation requires the stimulation current to penetrate the resistance of the skin and the electrode, about 1000 ohms. CEFAR REHAB X2 can penetrate this resistance and maintain a current of up to 99.5 mA. With a change in load from 100 to 1000 ohms, the stimulation current changes less than 10% from the set value.

The stimulator operates on either two non-rechargeable 1.5 V AA batteries or two rechargeable 1.2 V AA batteries, recharged in a separate charger.

### CEFAR REHAB X2

Number of channels .....	2 (non-independent)
Constant current .....	Up to a resistance of 1000 ohm (increased load can reduce the maximum current)
Stimulation current/channel .....	0–99,5 mA (charge maximum: 40µC)
Waveform .....	Symmetrical biphasic pulse, 100% compensated
Number of preset programs .....	27
Number of custom programs .....	3
Stimulation forms .....	Conventional Burst Modulated frequency/pulse duration Mixed frequency Alternated modulated pulse duration (CEFAR Flow TENS) Intermittent
Max pulse duration .....	400 µs
Max frequency .....	120 Hz
Timer .....	1 to 99 min./Off
Environment for storage, use and shipping .....	Temperature 10° C–40° C Air humidity 30%–75% Air pressure 700 hPa–1060 hPa
Power source .....	2 x 1.5 V AA non-rechargeable or 2 x 1.2 V AA rechargeable batteries
Current consumption for one channel, 80 Hz, 30 mA .....	150 mA
I r.m.s. max/channel .....	31 mA
Size .....	120 x 50 x 30 mm
Weight .....	ca. 180 g

# KEY TO THE SYMBOLS

 Reference number, part number

 Follow instruction for uses

 Type BF applied parts

 Council Directive 2012/19/EU concerning Waste Electrical and Electronic Equipment (WEEE) requires not to dispose of WEEE as municipal waste. Contact your local distributor for information of the unit and accessories

 Keep the device dry

 Minimum and maximum temperature indications to respect

 This equipment complies with all requirements of the Medical Device Directive (93/42/EEC)  
2797

 Manufacturer name and address, manufacturing year

 Power/Pause

## 14. ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The CEFAR REHAB X2 Electrotherapy System is intended for use in the electromagnetic environment defined below. The customer or the user of the CEFAR REHAB X2 Electrotherapy System should assure that it is used in such an environment		
Emissions Tests	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Group 1	The CEFAR REHAB X2 Electrotherapy System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	CEFAR REHAB X2 Electrotherapy System is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes
Harmonic emissions IEC 61000-3-2	Not Applicable - Battery powered	
Voltage fluctuations IEC 61000-3-3	Not Applicable - Battery powered	

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
The CEFAR REHAB X2 Electrotherapy System is intended for use in the electromagnetic environment specified below. The customer or the user of the CEFAR REHAB X2 Electrotherapy System should assure that it is used in such an environment			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Risk assessment on the CEFAR REHAB X2 Electrotherapy System indicates the compliance levels claimed are acceptable when ESD-precautionary measures are taken.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	Not Applicable - Battery powered  Not Applicable - signal lines less then 3 meters	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+1kV differential mode (line to line) +2kV common mode (line to ground)	Not Applicable - Battery powered	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines  IEC 61000-4-11	<5% $U_T$ (>95% dip in $U_T$ ) for 0,5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	Not Applicable - Battery powered	Mains power quality should be that of a typical commercial or hospital environment. If the user of the CEFAR REHAB X2 Electrotherapy System requires continued operation during power mains interruptions, it is recommended that the CEFAR REHAB X2 Electrotherapy System be powered from an uninterrupted power supply.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: $U_T$ is the AC mains voltage prior to application of the test level.			

### Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The CEFAR REHAB X2 Electrotherapy System is intended for use in the electromagnetic environment defined below. The customer or the user of the CEFAR REHAB X2 Electrotherapy System should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
<p>Conducted RF IEC 61000-4-6</p>	<p>3 Vrms 150 kHz to 80 MHz</p>	<p>[V] V, where V1 = 3V</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the CEFAR REHAB X2 Electrotherapy System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> $d = \frac{35 \sqrt{P}}{ V1 }$
<p>Radiated RF IEC 61000-4-3</p>	<p>3 V/m 80 MHz to 2,5 GHz</p>	<p>[E1] V/m, where E1 = 3V/m</p>	<p><math>d = \frac{35 \sqrt{P}}{ E1 }</math>    80 MHz to 800 MHz</p> <p><math>d = \frac{7 \sqrt{P}}{ E1 }</math>    800 MHz to 2,5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: center;">  </div>

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CEFAR REHAB X2 Electrotherapy System is used exceeds the applicable RF compliance level above, the CEFAR REHAB X2 Electrotherapy System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the CEFAR REHAB X2 Electrotherapy System.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.

Recommended separation distances between portable and mobile RF communications equipment and the CEFAR REHAB X2 Electrotherapy System			
<p>The CEFAR REHAB X2 Electrotherapy System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the CEFAR REHAB X2 Electrotherapy System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the CEFAR REHAB X2 Electrotherapy System as recommended below, according to the maximum output power of the communications equipment.</p>			
Rated maximum output power of transmitter P (W)	Separation distance according to frequency of transmitter d (m)		
	150 kHz to 80 MHz $d = \frac{3,5}{[V]} \sqrt{P}$ (where V1 = 3V)	80 MHz to 800 MHz $d = \frac{3,5}{[E]} \sqrt{P}$ (where E1 = 3V/m)	800 MHz to 2,5 GHz $d = \frac{7}{[E]} \sqrt{P}$ (where E1 = 3V/m)
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

**FCC REQUIREMENTS**

Part 15 of the FCC Requirements	
<p>This device complies with Part 15 of the FCC Rules. Operation is subject to the following 2 conditions:</p>	<ul style="list-style-type: none"> <li>• This device may not cause harmful interference</li> <li>• This device must accept any interference received, including the interference that may cause undesired operation.</li> </ul>
FCC ID	T9J-RN42
Contains Transmitter Module IC	6514A-RN42

# CLINICAL GUIDE

## INTRODUCTION


This clinical guide gives tips and practical advice on how to use the CEFAR REHAB X2 to treat common orthopedic and neurological indications.

The guide is organized by indication area. For each indication, the guide gives the goal of treatment and suggestions for where to place the electrodes. NMES is an active training method, and the pictures show how NMES can be integrated with training. For each indication, appropriate programs are suggested, but custom programs can also be created for patient-specific treatment.

A program with a short pulse duration, 200  $\mu$ s, is ideal for treating smaller muscle groups, while a longer pulse duration, 350–400  $\mu$ s, is better for larger muscle groups. The length of the work time and the rest time varies from program to program. For weaker muscles, we recommend programs with longer rest times, to give the muscle time to rest between contractions.

Just as with other training methods, NMES treatment should be adjusted to the patient's ability. Depending on the muscular status and where the patient is in the rehabilitation process, treatment can last from 5 to 60 minutes and be repeated from three times a week to twice a day. Remember that the patient may develop sore muscles after NMES treatment.

Suitable electrode sizes are suggested for each indication. We use the electrodes sizes shown below. The electrodes we recommend are shadowed.

- 1  $\varnothing$  32 mm
- 2  $\varnothing$  50 mm
- 3 40 x 60 mm
- 4 50 x 100 mm
- 5 80 x 130 mm
- 6 50 x 50 mm
- 7 50 x 90 mm

## UPPER ARM – REDUCED FUNCTION IN BRACHIAL BICEPS

### Examples of indications

1. Reduced function in the brachial biceps, for example after a fracture or spinal injury.
2. Reduced supination of the elbow.

### Goal



1. Increased strength and endurance in the brachial biceps.
2. Increased supination of the elbow.

### Placement of electrodes

Place one pair of electrodes over the brachial biceps.

### Suggested programs

CEFAR REHAB X2: 9, 11, 13, 19

## UPPER ARM – REDUCED FUNCTION IN BRACHIAL TRICEPS

### Examples of indications

1. Reduced function in the brachial triceps, for example after a stroke.
2. Elbow fracture.

### Goal




1. Increased strength in the brachial triceps and increased control of the elbow.
2. Increased strength in the brachial triceps.

### Placement of electrodes

Place one pair of electrodes over the brachial triceps.

### Suggested programs

CEFAR REHAB X2: 9, 11, 13, 19

**UPPER ARM – ALTERNATING STIMULATION**

**Examples of indications**


Extension/flexion defect in the elbow.

**Goal**

- 1. Increased mobility of the extension/flexion of the elbow.
- 2. Reduced spasticity.


**Placement of electrodes**

Place one pair of electrodes over the brachial triceps and one pair over the brachial biceps.

**Suggested programs**

CEFAR REHAB X2: 17, 18

**LOWER ARM – REDUCED STRENGTH/ HYPOTROPHY IN THE HAND EXTENSORS**

**Examples of indications**

- 1. Reduced wrist extension and/or flexion spasticity for example after a stroke.
- 2. Reduced wrist extension after fracture.

**Goal**


- 1. Facilitation and reduced spasticity.
- 2. Increased mobility of the wrist and increased strength in the hand extensor muscles.

**Placement of electrodes**

Place one pair of electrodes over the hand extensor muscles.

**Suggested programs**

CEFAR REHAB X2: 9, 11, 13, 15

### LOWER ARM – ALTERNATING STIMULATION

#### Examples of indications



Reduced extension and flexion of the wrist.

#### Goal

1. Increased extension and flexion of the wrist.
2. Reduced spasticity


#### Placement of electrodes

Place one pair of electrodes over the hand extensor muscles and one over the hand flexor muscles.

#### Suggested programs

CEFAR REHAB X2: 17

### HAND – THUMB EXTENSION

#### Examples of indications

1. Reduced extension of the thumb.
2. Flexor spasticity in the thumb.

#### Goal

1. Increased extension of the thumb.
2. Reduced spasticity in the thumb.


#### Placement of electrodes

Place a small oval electrode over the extensor pollicis brevis and a small round one proximal to the oval one.

#### Suggested programs

CEFAR REHAB X2: 9, 11, 13, 15

### HAND – REDUCED OPPOSITION/GRIP ABILITY

#### Examples of indications


Reduced opposition/grip ability.

#### Goal

Increased opposition and improved grip.

#### Placement of electrodes

Place a small electrode over the thenar muscles and another small one over the hypothenar muscles.

#### Suggested programs


CEFAR REHAB X2: 9, 11, 13, 15

**TRUNK, THORACIC SPINE – HYPOMOBILITY**

**Examples of indications**

- 1. Thoracic hypomobility.
- 2. Increased thoracic kyphosis.

**Goal**


- 1. Mobility training of the thoracic spine.
- 2. Improved posture.

**Placement of electrodes**

Place two pairs of electrodes along the spinal erector muscles paravertebrally and over the rhomboid muscles at the height of the hypomobile area.

**Suggested programs**

CEFAR REHAB X2: 10, 12, 14



**TRUNK, LOWER BACK – AWARENESS TRAINING OF BACK EXTENSORS**

**Examples of indications**

Pain/insufficiency in lower back.

**Goal**

- 1. Increased awareness of the trunk muscles.
- 2. Increased postural control.
- 3. Increased stability.


**Placement of electrodes**

Place two pairs of electrodes paravertebrally in the lumbar region, along the spinal erector muscles on both sides.

**Suggested programs**

CEFAR REHAB X2: 10, 12, 14

### KNEE – INSTABILITY

**Examples of indications**

1. Postoperative after ACL/knee operation.
2. Extension defect.
3. Post-stroke treatment.

**Goal**



1. Increased control of the quadriceps.
2. Increased mobility and decreased pain in the knee.

**Placement of electrodes**

Place one pair of electrodes over the quadriceps.

**Suggested programs**

CEFAR REHAB X2: 12, 14, 16

### KNEE – MUSCULAR IMBALANCE

**Examples of indications**

Muscular imbalance in the medial vastus in relation to the lateral vastus.

**Goal**



Increased strength in the medial vastus.

**Placement of electrodes**

Place one pair of electrodes over the medial vastus.

**Suggested programs**

CEFAR REHAB X2: 10, 12, 14

## FOOT – HALLUX VALGUS

### Examples of indications

Hallux valgus.

### Goal


Increased abduction of the big toe.

### Placement of electrodes

Place one pair of electrodes over the abductor hallucis.

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© 06/2020 DJO - 2948 - Rev C