# FRONT

## Pedal Exerciser with Digital Display Use Instructions

P554-00

### SAFETY INSTRUCTIONS:

- Never put your full body weight on the pedal exerciser.
- Check for loose, damaged, worn or missing parts during assembly and before. Do not use if any of these conditions exist.
- Regularly check for tight fit of the screws and nuts.
- Consult your doctor before starting this or any exercise program.

**TO UNFOLD:** Minding your fingers, push the leg away from body of the pedal exerciser until you hear the button click into place. Repeat with other leg.

**TO FOLD:** Minding your fingers, push the button and fold leg towards body of pedal exerciser. Repeat with other leg.

**SETTING THE RESISTANCE:** Screw resistance knob counterclockwise to loosen, clockwise to tighten.

### INSTRUCTIONS FOR USE (LEG EXERCISE):

1. Sit down on a stable object (e.g. a chair or the side of the bed).
2. Unfold the pedal exerciser and set in front of you on a firm, level and non-slip surface.
3. Set the pedal exerciser to a light resistance (screw resistance knob counterclockwise).
4. Begin pedaling slowly.
5. Pedal for short periods increasing speed, time intervals and resistance as you gain strength and endurance.
6. Keep good posture while pedaling (back and neck straight).
7. Breathe evenly and steadily in and out while pedaling.

**INSTRUCTIONS FOR USE (ARM EXERCISE):** As above but with pedal exerciser on a table.

### USE OF DIGITAL DISPLAY:

1. Activate display by pressing the red start button or by moving the pedal.
2. Scroll through the following display functions by pressing the red button. An arrow next to each function will indicate which is being displayed.
  - a. Time: Current work out time is displayed.
  - b. Count: Number of pedal revolutions is shown.
  - c. Calories: Approximate number of calories burned is displayed.
  - d. Scan: In this mode the display changes automatically between the different display functions giving you a successive overview of all workout results.

The display switches off automatically approximately 4 minutes after stopping.

### CLEANING AND MAINTENANCE:

The body and pedals can be cleaned with a damp cloth. Do not use harsh cleaners and do not let water enter the equipment.

Replace the AAA battery when the display cannot be read clearly anymore.

**STORAGE:** Store pedal exerciser in a dry place. Cover to keep dust free if stored for long period of time. Fold up to save space.

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All Rights Reserved