

CAREX[®]

ITEM# A33300 0000

Trio Rolling Walker



 **Please read all instructions before assembly or use.**

Questions, problems or missing parts?

Before returning, contact Compass Health Brands customer service:

800-526-8051 | CAREX.COM

WARNINGS

INDICATIONS FOR USE

This device is intended to provide ambulatory assistance for end users with moderate stability or balance issues.

SAVE THESE INSTRUCTIONS. This sheet must accompany device with each use.

CAUTION! Read all instructions and warnings prior to assembly and use.

- If you do not understand these instructions, contact a healthcare professional before attempting to use.
- Inspect all parts for shipping damages before use. If there is shipping damage – **DO NOT USE**. Contact the manufacturer for further instructions.
- If there are any problems with your rollator, **DO NOT** attempt to repair it yourself. Contact the manufacturer for any needed parts or repair instructions.
- Tires and brakes are non-durable components which will require maintenance through the life of the rollator.
- Brakes should be adjusted so that the wheels do not move while the brake is applied – if brakes are unable to be adjusted correctly, **DO NOT** continue to use the rollator, seek repair or replacement.
- Tires should be inspected regularly for damage and wear and replaced when the tires are excessively worn. Tires may wear faster if the device is used outdoors on rough surfaces.
- Compass Health Brands assumes no responsibilities for any damage or injury caused by improper installation, assembly, or use of this product.


INSTRUCTIONS FOR USE AND WARNINGS

- This rollator has a weight capacity of 300 lbs. (136 kg). **ALWAYS** observe the weight capacity on the labeling of your rollator.
- **ALWAYS INSPECT ROLLATOR PRIOR TO USE.** Properly inspecting and maintaining your rollator can reduce the risk of injury. Make sure that all parts are secure and are in good working order. **ALWAYS** check the fasteners, wheels and brakes for tightness prior to use. Make sure that parking brakes are working, and the rollator wheels cannot turn with the brakes on. If rollator is not in good working order or fasteners cannot be tightened, **DISCONTINUE USE OF THE ROLLATOR.**
- All wheels and handles **MUST** be adjusted to the same height for safe use. **DO NOT** use the rollator if wheels or handles are out of alignment, or if hardware is not secure. As applicable to your rollator model, all hand knobs must be firmly tightened and height adjustment buttons must protrude fully through the height adjustment holes with an audible “click”.
- **DO NOT** use this device if you are unable to evenly balance your weight between the two handles. **DO NOT** use this device if you have a history of dizziness, balance problems, or a need for a device that can support heavy leaning. Consult a healthcare provider to determine if this device is right for you.
- **ALWAYS** wear protective footwear while using the rollator. **ALWAYS** keep shoelaces, loose clothing and/or accessories clear of the wheels and handle bars during use.
- **ALWAYS** maintain good posture while using the rollator. When walking with the rollator, keep your feet in line with the rear wheels, and **DO NOT** let the rollator get too far in front of you.
- **DO NOT** use the rollator to pull yourself up from a chair or other seated position. Use the arms of the chair to support your weight as you stand.
- **ALWAYS** have all wheels in contact with the floor/ground while using this device. Failure to maintain wheel contact could result in improper balance of the rollator and cause injury or damage.
- **ALWAYS** make sure the rollator is fully opened and the folding hinge is locked in the down position before use.
- **DO NOT** hang anything from the rollator handles or frame. Items may only be carried in the provided pouch/basket.
- **DO NOT** climb, kneel, or stand on the rollator. This is not a toy. **DO NOT** let children play on rollator.
- **BE AWARE** of your surroundings when operating the rollator. Look for hazards and avoid them. Avoid streets and surfaces with water, sand, gravel, dirt, leaves, and other debris. Wet, slick, uneven and/or rough surfaces may impair traction and contribute to possible accidents. Loose cords and unsecured rugs may move suddenly and cause a loss of balance. **DO NOT** use a rollator in mud, ice, or puddles.
- **NEVER** use near steps, sloped driveways/hills with a steep incline (never greater than 15 degrees). **DO NOT** walk backwards down a slope.
- **DO NOT** use stairs or escalators while using the rollator. **DO NOT** use near swimming pool areas.
- **ALWAYS** take additional care and move slowly when moving from carpeted to hard surface floors or cracks/separations in walkway surfaces; avoid sharp bumps, drainage grates, and sudden surface changes. Rollator may abruptly stop if a wheel becomes wedged.
- **DO NOT** use the hand brake as a primary stopping tool. It may be used if needed but is primarily designed to maintain a stopped stance.
- **DO NOT** use while under the influence of mind-altering drugs, alcohol or while suffering from dizziness.
- **DO NOT** store in direct sunlight or outdoors. Exposure to extreme temperatures (heat, cold, moisture, etc.) may make surfaces (brake handles, etc.) hot to the touch, weak or brittle. Ensure brake handles **DO NOT** twist and skin contact surfaces are cool to the touch before use.

**FAILURE TO USE COMMON SENSE AND HEED THE ABOVE WARNINGS
FURTHER INCREASES RISK OF SERIOUS INJURY. USE AT YOUR OWN RISK AND WITH
APPROPRIATE AND SERIOUS ATTENTION TO SAFE OPERATION. USE CAUTION.**

1 Check Parts

- Unpack and check that you have all the parts and attachments shown.
- To avoid risk of suffocation to children or animals, dispose of the plastic bags correctly.



Description	Qty
Bolts	2
Washers	2
Hand Knobs	2
Threaded Hand Knobs	2
Frame	1
Handles	2

NOTE: Actual items may vary from what is shown.



Check each part carefully for shipping damage. **DO NOT** use the trio rolling walker if there are signs of damage and/or parts missing.



WARNING: To avoid serious injury, the trio rolling walker must be properly assembled. Ensure all bolts are fixed in place and hand screws have been firmly tightened before each use.

Product Assembly


2 Open Frame

- Open your Trio Rolling Walker. Push the folding hinge until straight and locked to secure the unit in the proper position.



3 Pull Out Rear Wheels

- Pull out rear wheels until you hear a "click".



4 Secure Rear Wheels

- Insert threaded hand knob and hand-tighten. Repeat step on the opposite side.



5 Assemble Handles

- Insert handles into frame.




6 Adjust Handle Height

- Adjust handle height as shown.
- Line up holes and secure with hex-head screw, washer and hand knob. Repeat with other handle.



7 Braking System

- Attach brake cable clamps.



Pull up to apply brakes



Release brakes



Push down to lock brakes

8 Brake Adjustment

- Brake cables may need adjusting.

To adjust brakes:

- Find nut near bottom of brake line
- Hold lower nut with needle nose pliers
- Adjust upper nut with 5/16" or 8 mm wrench
- Turn counter clockwise to loosen
- Turn clockwise to tighten



Cleaning and Maintenance

1. Clean the trio rolling walker using mild soap and water. Wipe dry.
 2. Your trio rolling walker should be checked periodically to ensure the wheel locks are functioning properly and that all hardware is secure.
 3. Check that the front castor wheel swivels freely. **DO NOT** use the trio rolling walker if it does not.
-

CAREX[®]

Questions, problems or missing parts?

Before returning, contact Compass Health Brands customer service:

800-526-8051 | CAREX.COM

Retain this manual for future use.

COMPASSHEALTH

Manufactured for:
Compass Health Brands Corp.
Middleburg Heights, OH 44130
www.compasshealthbrands.com
Phone 800.376.7263
Made in China
42-A33300_0000_05 © 2022