- Correctly re-apply the device to the patient's leg and ensure the risk of deep vein thrombosis and pulmonary embolism.

#### GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device. Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, excessive pressure, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

#### FITTING INSTRUCTIONS

The device is a truly custom brace adapted to the patient's anatomy, activity level, and desired configuration. Therefore, images within the Instructions for Use may appear different from the final product. Depending on the selected device option, please refer to either the SmartDosing® option (**Fig. 1**) or the Standard Dynamic Force System™ (DFS) option (**Fig. 2**) overview figure for locating components mentioned in the following instructions.

#### Device Application

- Open the Buckle (1A) and/or unfasten straps. Ask patient to sit down and extend leg while fitting the device. Place the device on the patient's leg with the Hinge (1B, 2A) on the affected side of the knee.
- Ensure proper alignment of the device on the leg.
- Height positioning: Align the center of the Hinge slightly above the middle of the patella.
- Side positioning: Center of the Hinge should be at midline of leg.
- Wrap the Upper Calf Strap (1C, 2B) around the back the leg above the calf muscle, thread it through the D-ring, and secure (**Fig. 3**).
- Hold the Thigh Shell (1D, 2C) and bring the Thigh Strap (1E, 2D) around leg. Thread the Thigh Strap through the slotted thigh strap attachment and secure (**Fig. 4**).

#### SmartDosing option (Fig. 1)

- Place button of the Buckle into appropriate keyhole and close carefully using the palm of the hand to snap the Buckle closed (**Fig. 5**).
- With the patient's knee fully extended, adjust the DFS Strap (1F) length until it sits firmly against the leg, and then fold it into Alligator Clip (1G). At this point, patient should not experience any tension or unloading.
- Ask patient to bend knee with the foot flat on the floor. Turn the SmartDosing Dial (1H) clockwise until the indicator is in the “5” position. Have the patient stand up and take a few steps to verify correct positioning of the device and tightness of the straps.
- Determine optimal DFS Strap tension based on patient's pain relief feedback.
- If patient requires more or less tension with indicator at the “5” position, adjust the length of the DFS Strap accordingly.
- Aim for final SmartDosing Dial setting at the “5” position since this will give patient the ability to adjust dosing during daily life activities.

#### Standard DFS option (Fig. 2)

- Thread the Lower Calf Strap (2E) through the D-ring and secure (**Fig. 7**).
- Bend the patient's braced knee to an approximate angle of 80 degrees. While holding the Thigh Shell, to prevent the device from moving, thread the DFS Strap (2F) through D-ring in front of the Thigh Shell and tighten (**Fig. 8**). The DFS Strap should be worn snug for maximum pain relief.
- Determine optimal strap tension based on the patient's pain relief feedback.
- If patient requires more or less tension, adjust the length of the DFS Strap accordingly.

#### Device Adjustments

#### Varus and valgus angle adjustments

The device's Adjustable Dynamic Joint™ (ADJ) hinge (**Fig. 9**) can be used to adjust varus and valgus angles, as well as to create an offloading force to either the medial or lateral compartment. For one ADJ mechanism at a time, carry out the following instructions to adjust the angles:

- Remove the ADJ cover to reveal the adjusting mechanism.
- Turn the Locking Screw (A) counterclockwise approximately two full turns to allow for angle adjustments (do not remove completely).
- Turn the Adjusting Screw (B) clockwise to increase hinge angle or counterclockwise to decrease hinge angle.
- Half a turn of the Adjusting Screw represents approximately one degree adjustment. Therefore, one complete turn of the screw equals approximately two degrees of adjustment.
  - Note:** Over-adjustment may result in brace rotation and/or excessive pressure. Discrepancy between proximal and distal correction may result in altered fit of the device.
- When the desired angle has been achieved, retighten the Locking Screw.
  - Note:** To prevent Locking Screw from loosening, a thread locking compound may be applied.
  - Note:** Do not overtighten the Locking Screw.
- Repeat the adjustment procedure for the other ADJ mechanism.
- When complete, replace the ADJ covers.

7. Correctly re-apply the device to the patient's leg and ensure the fit. Ensure the DFS Strap is applied with sufficient tension to augment the unloading effects of the ADJ hinge.

#### Range of motion (ROM) adjustments

All hinge options are shipped with a pre-selected ROM stop in place.

#### Poly-axial hinge

If the ROM needs to be adjusted for the device with a mounted Poly-axial hinge, please contact Ossur Customer Service or Technical Department for assistance.

#### Device Removal

1. Ask patient to sit down with the leg extended.

#### SmartDosing option (Fig. 1)

- Turn SmartDosing Dial counterclockwise until the indicator is at the “0” position to release tension on the DFS Strap.
- Bend the patient's knee to 90° and unfasten the Upper Calf Strap, Thigh Strap, and the Buckle. Pull the Buckle button out of Keyhole.

#### Standard DFS option (Fig. 2)

- Remove the DFS Strap from the D-ring.
- Unfasten the Thigh, Upper, and Lower Calf Strap from the outer side of the thigh.

#### Accessories and Replacement Parts

Please refer to the Ossur catalog for a list of available replacement parts or accessories.

#### USAGE

#### Cleaning and care

Washing the device with the soft goods detached allows for more thorough cleaning.

#### Washing Instructions

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

**Note:** Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

**Note:** Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

#### Hinge

Remove foreign materials (e. g., dirt or grass) and clean using fresh water.

#### DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

#### LIABILITY

Ossur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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