ENGLISH

SYMBOLS

MD Medical Device

Caution: This device contains natural rubber **latex** which may cause allergic reactions.

INTENDED USE

The device is intended for uncompartmental unloading of the knee

Indications for use

For osteoarthritis and other knee conditions that may benefit from unloading and/or compression, such as:

- Medial or lateral knee pain
- Early or mild uni-compartmental osteoarthritis of the knee
- Symptomatic, uni-compartmental /meniscus degeneration

Contraindications

Must not be used by individuals for whom compression is contraindicated.

Warnings and Cautions:

- The device should not be used for extended periods of time and removed during long periods of rest, e.g., sleep.
- Patients with already compromised blood circulation in their extremities, like those with diabetes or peripheral vascular disease, need to be extra careful and are advised to consult with a healthcare professional before using the device. If you experience pain, swelling, sensation changes, or if your extremity shows signs of insufficient blood flow (turns blue, white, or cold) while using this device, discontinue the use immediately and consult with a healthcare professional.

Device Application

- First loosen the two ends of the DFS straps without pulling them out of the slits.
- For easy donning, fold the top part of the support outward so that the FoldLock edge becomes visible. Grasp the support with two hands on the side and pull the support over the knee with the hinge on the affected side of the knee (**Fig. 1**).
- Pull evenly over the knee joint until the knee cap is free from pressure and centered into the patella opening. Then fold the edge of the support back into the original position (**Fig. 2**).
- Bend the knee to an 80° angle and tighten the DFS straps until pain relief is sufficient (**Fig. 3**). With the knee in full extension, verify correct position and strap tension. Readjust, if needed (**Fig. 4**).

Device removal
Fold the FoldLock edge of the support outward so that the edge becomes visible. Undo the DFS straps. Attach the hook and loop straps appropriately to avoid damaging the knit. Grip the support at the level of the calf and pull it off in a downward direction.

Cleaning and care

Wash regularly to retain the compression of the textiles.

- Machine-wash with a mild detergent at cold temperature (30°C)
- Air dry

Note: Secure the hook and loop fasteners before washing to avoid damaging the textiles.

Note: Do not tumble dry, iron, bleach, or wash with fabric softener.
Note: If the device is used in salt water or chlorinated water, rinse with fresh water and air dry.

DISPOSAL
The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

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LIABILITY
Osour does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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