














Instructions for Use

REBOUND[®] DUAL ST

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	10
FR	Notice d'utilisation	16
ES	Instrucciones para el uso	21
IT	Istruzioni per l'uso	26
NO	Bruksanvisning	31
DA	Brugsanvisning	36
SV	Bruksanvisning	41
EL	Οδηγίες Χρήσης	46
FI	Käyttöohjeet	52
NL	Gebruiksaanwijzing	57
PT	Instruções de Utilização	62
PL	Instrukcja użytkownika	67
CS	Návod k použití	73
TR	Kullanım Talimatları	78
RU	Инструкция по использованию	83
JA	取扱説明書	89
ZH	中文说明书	94
KO	사용 설명서	99

Medical Device

PRODUCT DESCRIPTION

This document covers the following product versions:

- Rebound Dual ST
- Rebound Dual ST Custom
- Rebound Dual ST Custom with Acculign

List of Parts (Figure 1):

1. Lower Shin Strap #1
2. Upper Thigh Strap #2
3. Lower Thigh Strap #3
4. Upper Shin Strap #4
5. Calf Strap #5
6. Upper Frame
7. Lower Frame
8. D ring
9. Hinge, Condyle pad
10. PCL Slot

List of Parts (Figure 14)

11. Acculign™ Hinge Arm
12. Acculign™ Hinge Arm Screw

INDICATIONS FOR USE

The brace is indicated for combined and non-combined instabilities of the:

- Anterior cruciate ligament (ACL)
- Medial collateral ligament (MCL)
- Lateral collateral ligament (LCL)
- Posterior cruciate ligament (PCL)

Acculign version only: Conditions requiring unicompartmental load control.

CONTRAINDICATIONS

None known.

INTENDED USE

Target group includes individuals requiring extra knee joint support, individuals with ACL, MCL, LCL, PCL and combined instabilities, participants in low to high impact activities and sliding sports, and everyday athletes.

GENERAL SAFETY INSTRUCTIONS

This knee brace is a supportive knee brace only. It is not intended or guaranteed to prevent knee injury.

The purchaser agrees to the following:

1. Össur accepts no responsibility for any injury which might occur while using Össur knee bracing products.

2. With use of this type of device, it is possible to receive various types of skin irritation, including rashes and pressure points. If you experience any unusual reactions while using this product, consult your health care provider immediately.
3. If you have any questions as to how this knee brace applies to your own specific knee condition, consult your health care provider.
4. If brace causes discomfort, irritation, or is damaged, discontinue use immediately and contact your health care provider.
5. This brace is single patient use only. If you do not agree to the above conditions, contact your brace provider before using this product.

FITTING INSTRUCTIONS

Brace Application

1. Unfasten straps 1#,2#,3# and 4#. Loosen strap #5 (**Figure 2**).
2. Bend your knee to between 30 - 45 degrees, step through the brace, and then slide it up your leg. Align the Össur hinge alignment guide (middle of the hinge) with top of the knee cap. Make sure the condyle pads are in full contact with the knee (**Figure 3**).
3. Insert lower shin strap #1 through the D ring and fasten tightly (**Figure 4**).
4. Insert upper thigh strap #2 through the D ring and fasten (**Figure 5**).
5. Then fasten lower thigh strap #3. Tightening the lower thigh strap #3 will pull the hinges towards the back of the knee (**Figure 6**). **note:** Do not overtighten the two top straps as this can make the brace move.
6. Insert upper shin strap #4 through the D ring and tighten until the hinges are aligned with the midline of the knee. Readjust strap #3 if necessary (**Figure 7**).
7. Tightly fasten calf strap #5, while making sure it fits comfortably. After fitting, the calf strap will only need to be readjusted to improve brace comfort/ suspension (**Figure 8**).
8. When standing, the brace should fit the leg with no larger gaps. The hinge should be aligned with the middle of the knee cap. Seen from the side, the hinges are aligned with the knee midline (**Figure 9**).

TROUBLESHOOTING

Proper anterior/posterior placement will help prevent migration and ensure proper brace suspension.

- If the brace is too far in the front, loosen the upper shin strap #4, and then tighten the calf strap #5 and lower thigh strap #3.
- If the frame is too posterior, loosen the calf #5 and lower thigh #3 straps and tighten the upper shin strap #4.
- Pay attention to ensure the hinge is positioned correctly. Incorrect positioning may result in the brace slipping down.
- Incorrect positioning or fastening the straps in the wrong order may result in slippage of the brace.

BRACE ADJUSTMENTS

Adjusting Pad and Strap Length:

For optimal fit and comfort, ensure the strap pads are centered on the leg when tightening straps. If a strap pad is too long, remove the pad from the strap and use scissors to shorten the pad. To adjust the length of a strap, detach the alligator clip at the end of the strap. Adjust the length as necessary and the alligator tab to the strap.

Widening or narrowing fit at knee joint:

The brace can be adjusted to individual knee widths. Adjust to knee width by removing the medium condyle pads and replacing them with the thick or thin version. Any combination of pads can be used to achieve a more comfortable fit.

Adjusting the Brace Contouring:

The Rebound DUAL ST brace can be adjusted to individual anatomies. Use bending irons to cold form the aircraft grade aluminum frame to fit different thigh and calf shapes.

Replacing and Reconfiguring D rings:

D rings are field-serviceable. D rings can be moved from the interior to the exterior or replaced if broken.

To move a D ring from the interior to the exterior:

1. Remove strap from D ring.
2. Bend D ring mounted on the brace interior and push it inwards to release from the keyhole lock. Remove from the rivet (**Figure 10**).
3. Attach D ring to the brace exterior by pushing down onto rivet through the keyhole lock. Slide and click D ring into place. Place keyhole lock with indented side facing up (**Figure 11**).
4. Replace strap. If needed, repeat steps for other D rings.

Adjusting the Varus / Valgus angle of the Brace (Acculign version only):

The brace can be adjusted for individual varus or valgus angles and to create an offloading force to either the medial or lateral compartment.

- To adjust the angle, loosen the screws on the desired side of the brace and slide the upper frame to the desired position with a steady force.
- Adjust the medial upper frame to offload the medial compartment or adjust the lateral upper frame to offload the lateral compartment. The slide bar has a lettered scale and may be moved to increase or decrease the amount of varus or valgus force on the leg. (**Figure 14**). The lettered scale does not relate to degrees of offloading and should be used as a guide.
- Retighten the screws until snug (recommended torque value of 3 in-lb).
- All adjustments should be determined by a healthcare provider, in accordance with the patient's comfort level. Each time the brace is adjusted, the healthcare provider should note the letter for a record of progression.

Adjusting the Height of the Brace (Acculign version only):

- The brace can be adjusted for individual height requirements.
- The brace height can be adjusted by loosening the screws on both sides of the upper frame and pulling the frame up with a steady force. Ensure the brace is properly aligned by checking that both hinge arms are located on the same position on the medial and lateral hinge arm.
- Retighten the screws until snug.

Range of Motion adjustment (ROM) (optional accessory, not included):

The brace has 0° extension stops attached to the hinges. If needed, range of motion can be adjusted using the extension and flexion stop kits.

Extension and flexion stop kits are available separately for order.

To install a kit:

1. Remove the condyle pads from the medial and lateral hinges. Remove retaining screws.
2. Using a short Philips head screwdriver, slide the 0° stop out of the hinge on both sides.
3. Re-apply the desired extension stop and gently secure it with the retaining screw.
4. Repeat steps to change the flexion stop.

NOTE: Do not overtighten retaining screw.

PCL Strap Attachment (optional accessory, not included):

The PCL kit is designed to counteract posterior subluxation of the tibia. Trim and adjust the PCL strap for proper fit.

Attaching PCL Strap – Acculign version

1. Note height position on medial and lateral hinge arms.
2. Remove hook-and-loop fastener to access connection area.
3. Remove Acculign hinge screws.
4. Remove Acculign hinge d-Ring.
5. Position PCL Acculign hinge d-ring.
6. Position PCL Acculign hinge inserts with raised portion into recessed pocket.
7. Return frame to initial positions on hinge arms.
8. Tighten Acculign screws.
9. Reapply hook-and-loop fastener.
10. Attach PCL strap.

Attaching PCL Strap – other versions

1. Take the end of the PCL strap with the alligator clip and loop it through the PCL slot on the brace frame. Alligator clip will close on the inside of the strap. Trim strap if necessary (**Figure 12**).
2. Loop the other end of the PCL strap through the frame. Fasten strap and secure it onto the outside of the strap (**Figure 13**).
3. Place the adjustable pad in the middle of the strap. Trim the strap pad if needed.
4. Fully extend the leg. Tighten the strap, securing the pad.

PRODUCT CARE

Follow the steps below to keep your brace in a good condition. Liners, straps and pads can be removed for cleaning. Hand-wash using mild non-detergent soap. Rinse thoroughly, Air-dry.

If the brace moves, clean lower liners and pads with an alcohol wipe or damp cloth.

NOTES:

- Do not machine wash
- Do not tumble dry

Hinge Maintenance

Regularly check the hinge regularly for dirt, debris or similar contamination. Remove dirt, debris or similar contamination carefully using water or compress air.

NOTE: Do not disassemble the hinge.

REPLACEMENT PARTS

- Condyle Pads
- Liner Kits
- Strap Kits
- Color Pop Kits

OPTIONAL ACCESSORIES

- Range-of-motion (ROM) Kit
- PCL Kit

REPORT A SERIOUS INCIDENT

Important notice to users and/or patients established in Europe:

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability if the device is not maintained as instructed by the instructions for use, if the device is assembled with components from other manufacturers, or if the device is used outside of recommended use condition, application, or environment.

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