# COMPEX®

*Fit1.0*

*Fit3.0*

*SP2.0*

*SP4.0*

## **QUICK START GUIDE**

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*Beknopte handleiding*

*Guia rápido*

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## ***RULES FOR USE***

### **1. POSITION OF ELECTRODES**

For optimum results, use the recommended electrode positions.

A different electrode must be connected to each end of the cable.

N.B. It is occasionally the case that large dual-connection electrodes are used in such a way that a connection is not made.

It is important to select appropriately sized electrodes (small or large) and to position them correctly on the muscle group that you wish to stimulate in order to ensure effective treatment. It is therefore appropriate to always use the electrode size shown in the illustration. If necessary, try to find the best position possible by moving the electrode on the muscle, until you locate the point that produces the best contraction, or the greatest comfort is found.

For TENS type programmes, the general rule is to cover the painful area with electrodes.

N.B. Approximate positioning of electrodes makes the session less effective, but is not dangerous.

### **2. BODY POSITION DURING STIMULATION**

This position varies according to the muscle group that you wish to stimulate and the programme that you are using.

For programmes inducing powerful muscular contractions (tetanic contractions), it is recommended to stimulate the muscle isometrically. You should therefore hold the extremities of limbs being stimulated.

This position opposes maximum resistance to movement and avoids any shortening of the muscle during contraction, which would risk causing cramps. For example, when stimulating the quadriceps, you should be in a sitting position with ankles held to avoid extending the knees.

For the other types of programmes which do not induce powerful muscular contractions but only muscular twitches or tingling, position yourself as shown on the illustration, maintaining a comfortable position.

### **3. ADJUSTING STIMULATION INTENSITIES**

On a stimulated muscle, the number of fibres used depends on the stimulation intensity.

For programmes inducing powerful muscular contractions (tetanic contractions), you must use the maximum stimulation intensities (up to 999) while remaining within your tolerance threshold, in order to use the maximum number of fibres. The progress of a muscle will be greater the higher number of fibres that are worked. During a session it is therefore important to try to increase the intensity every 3 to 5 contractions.

For other types of programmes like recovery, massage, capillarisation or even muscle pain which only induce muscle twitches, you should gradually increase the stimulation intensities until clearly visible muscle twitches are obtained.

For TENS, Epicondylitis and Tendinitis programmes, you should gradually increase the intensities until you feel a tingling sensation under the electrodes.

### **4. PROGRESSION OF STIMULATION CYCLES**


It is recommended to start with the 1st cycle and progress to the next level when the cycle is complete, normally after 4 to 6 weeks of stimulation based on 3 sessions per week. It is also important to have reached significant stimulation intensities in sessions before going on to another cycle.

At the end of the cycle you can either start a new cycle or do maintenance training based on one session per week.

## START USING YOUR COMPEX

To help you get started with your Compex stimulator, connect to the website: [www.compex.com](http://www.compex.com)

### 1. DEVICE DESCRIPTION

**A** On/off button

**B** Sockets for the 4 stimulation cables

**C** Stimulation cables

**D** +/- buttons for the 4 stimulation channels

**E** I-button, which allows:


- The intensities to be increased on several channels at the same time
- The last 5 programmes used to be accessed

**F** Charger socket (slide the red cover to the right to reveal the charger connector)

**G** Battery compartment

**H** Socket for belt clip

## 2. CONNECTING THE CABLES

N.B. The MI-sensor cable (if available with the device) can be connected to any socket on the stimulator.

## 3. LANGUAGE, CONTRAST AND VOLUME SELECTION


- B** Press the +/- button of channel 1 to select the language to be used.
- C** Press the +/- button of channel 2 to adjust screen contrast.
- D** Press the +/- button of channel 3 to adjust the volume.
- E** Press the +/- button of channel 4 to adjust the backlighting.
- A** Press the on/off button to confirm and save your selections. Settings are immediately applied.

This options screen can then be displayed by switching the device off and by pressing and holding the on/off button for over 2 seconds.


## 4. CHOOSING A CATEGORY

N.B. The following screens are generic examples but they work in the same way regardless of the device that you have.

- A** Press the on/off button to switch the device off.
- B** Press the +/- button of channel 1 to select a category.
- E** Press the +/- button of channel 4 to confirm your selection.

## 5. SELECTING A PROGRAMME


- A** Press the on/off button to return to the previous screen.
- B** Press the +/- button of channel 1 to select a programme.
- E** Press the +/- button of channel 4 to confirm the selection and start the stimulation session.

N.B. After selecting a programme an additional screen may appear requesting the programme cycle to be selected. In either case, select the 1st cycle and then press the +/- button of channel 4 to confirm the selection and start the stimulation session.

## 6. LAUNCHING A PROGRAMME


N.B. If the MI-sensor cable is connected, the MI-scan test starts immediately after the programme is selected. At the end of the test, intensities must be increased to start the programme.

**A** Press the on/off button to switch the device into Pause mode.

**BCDE** The four channels flash, going from + to 000. The stimulation intensity must be increased to be able to start stimulation. For this, press the + buttons of the relevant channels until the desired setting is reached.

## 7. END OF A PROGRAMME

**A** To stop the stimulator, press the on/off button.

Please refer to the instructions for a more detailed description of the sequences presented in this quick start guide.

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