**DONJOY®**

## ULTRASLING® PRO AB and ER/IR

ULTRASLING® PRO AB Y ER/IR  
ULTRASLING® PRO AB UND ER/IR  
ULTRASLING® PRO AB ED ER/IR  
ULTRASLING® PRO AB ET ER/IR  
ULTRASLING® PRO AB EN ER/IR  
ULTRASLING® PRO AB E ER/IR  
ULTRASLING® PRO AB A ER/IR  
ULTRASLING® PRO AB A ER/IR  
ULTRASLING® PRO AB OG ER/IR  
ULTRASLING® PRO AB OCH ER/IR  
ULTRASLING® PRO AB JA ER/IR  
ULTRASLING® PRO AB OG ER/IR  
ULTRASLING® PRO AB 和 ER/IR  
ULTRASLING® PRO AB および ER/IR

ENGLISH .....	2
ESPAÑOL.....	4
DEUTSCH .....	6
ITALIANO .....	8
FRANÇAIS.....	10
NEDERLANDS .....	12
PORTUGUÊS.....	14
ČEŠTINA .....	16
SLOVENČINA.....	18
DANSK.....	20
SVENSKA.....	22
SUOMI .....	24
NORSK.....	26
简体中文 .....	28
日本語 .....	30

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
## ENGLISH

**BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.**

### INTENDED USER PROFILE:

The user should be able to:

Read and understand the directions, warnings and cautions.

**INTENDED USE/INDICATIONS:** Post-operative rotator cuff repairs, Bankart procedures, Capsular shifts, Glenohumeral dislocation/subluxation, Soft tissue strains/repairs.

**CONTRAINDICATIONS:** None.

**WARNINGS AND CAUTIONS:** If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, consult your medical professional immediately.

### SIZING AND CONFIGURATION:

**Arm Envelope:** To size the arm envelope, measure elbow crease to base of index finger.

ARM SLING SIZE	MEASUREMENT
S/M	10.5"-15" (26.7cm – 38.1cm)
L/XL	15+" (38.1+cm)

The arm envelope is shipped in L/XL configuration. To adjust arm envelope to the S/M size, unbuckle both quick release buckles. Fold in arm envelope and adhere the hook and loop tabs on inside of envelope.

**Right to Left configuration:** The UltraSling PRO is shipped in a right arm injury configuration. To adjust the shoulder ring to a left arm injury configuration, release the two quick release buckles at front of arm envelope and detach the cushion from the arm envelope. Apply cushion to left side of body by flipping cushion over so smaller end is pointing to back of left elbow. Wrap waist belt around waist and attach arm envelope to cushion. Slip the shoulder ring over the right shoulder. Place arm in sling and reattach the quick release buckles.

**Waist Belt:** The UltraSling PRO fits waist sizes 25" – 63" (63.5cm–160cm). Adjust waist belt and trim as necessary for proper fit.

**Optional pillow:** The UltraSling PRO can be worn without the pillow. Remove back of waist belt from pillow and adhere y-strap to both sides of the arm sling envelope at the elbow. For front of the waist belt, remove quick-release buckle from pillow and adhere y-strap to arm sling envelope. Trim the waist belt as needed.

### APPLICATION:

To ensure proper fit, if possible, product should be sized on patient pre-operatively.

1. Release quick-release buckle (i) on cushion and forearm strap (ii) on arm envelope. Unfasten the shoulder ring by releasing the quick release buckles (iii) at the front of the arm envelope. (Fig. A).
2. Align cushion on the injured side at waist/hip level with elbow fixed at 90°. Narrow end of cushion should be positioned to the back of forearm/elbow. Bring waist strap around back and attach to quick-release buckle at front of pillow. Adjust waist belt and trim as necessary for proper fit. (Fig. B).

**AB:** Place the large pillow under the affected area with the 45° angle facing up and wrap the waist belt around the waist, adhering the hook and loop to the pillow. Attach the sling to the cushion. (Fig B1). Adjust strap and trim as necessary for proper fit.

**ER/IR:** Align cushion on the injured side at waist/hip level with elbow fixed at 90°. Place cushion so line on top of cushion is parallel to the front of your body (if you stand at a table, the line should be parallel to the table). Wrap the waist belt around the waist,

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adhering the hook and loop to the pillow. Attach the sling to the cushion. (Fig B2). Adjust strap and trim as necessary for proper fit.

3. Place forearm into arm envelope as far back as possible. Secure thumb strap (iv) between the thumb and forefinger at front of sling. Secure forearm strap to sling and cushion. (Fig. C).
4. Wrap shoulder ring around healthy shoulder. Adjust shoulder ring circumference for proper fit. (Fig. D)  
**NOTE:** The hook and loop closure of the shoulder ring should be facing up and away from the patient.
5. Fasten the shoulder strap to the arm envelope by using the quick release buckles (iii). Adjust shoulder and back strap length for proper fit. (Fig. E).

#### **BRACE REMOVAL:**

To remove the UltraSling PRO, detach the quick release buckles, remove the forearm strap, thumb strap and waist belt. Slowly remove brace.

#### **CLEANING & STORAGE INSTRUCTIONS:**

	Hand wash in water (86°F/30°C), using mild soap. Air dry. <b>NOTE:</b> If not rinsed thoroughly, residual soap may cause irritation and deteriorate material.
	Do NOT iron.
	Do NOT tumble dry.
	Do NOT bleach.

**MATERIAL CONTENTS:** Nylon, Polypropylene, Polyurethane Foam, Cotton, Laminate Foam/Tricot, Polyacetal, Polyester, Acetal.

**WARRANTY:** DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale. To the extent the terms of this warranty are inconsistent with local regulations, the provisions of such local regulations will apply.

#### **RX ONLY**

**INTENDED FOR SINGLE PATIENT USE.**

**NOT MADE WITH NATURAL RUBBER LATEX.**

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A LICENSED HEALTH CARE PROFESSIONAL.

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