

ProCare® Quick-Fit™ Wrist & Forearm

Patient Application Instructions



APPLICATION INSTRUCTIONS:

BEFORE USING THE DEVICE, PLEASE READ THE **INSTRUCTIONS FOR USE** COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USE / INDICATIONS:

Treatment of wrist sprains, strains, and Carpal Tunnel Syndrome (CTS).

CONTRAINDICATIONS:

None

APPLICATION INFORMATION:

- 1 Loosen all straps and position hand in brace. Make sure the palmer stay is positioned comfortably below the palmer crease and aligned with forearm.
- 2 Position adjustable dorsal stay casing to ensure maximum support and proper fit. Stay casing may be moved proximal or distal for desired controlled range-of-motion.
- 3 Use single pull straps to secure body of brace. Straps should be snug but not constrict.
- 4 Secure hook-and-loop strap between thumb and forefinger.



WARNINGS AND PRECAUTIONS:

If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, consult your medical professional immediately. Never wear support over an open wound.

CLEANING INSTRUCTIONS:

Hand wash in warm water using mild soap, rinse thoroughly. AIR DRY. Note: If not rinsed thoroughly, residual soap may cause irritation and deteriorate material.

PATENT PENDING

FOR SINGLE PATIENT USE ONLY.

NOT MADE WITH NATURAL RUBBER LATEX.



Individual results may vary. Neither DJO, LLC nor any of its subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.

DJO, LLC | A DJO Global Company
T 800.336.6569 | D 760.727.1280 | F 800.936.6569
2900 Lake Vista Drive | Lewisville, TX 75067 | U.S.A.
DJOGlobal.com