	USERS MANUAL.....1
	BENUTZERHANDBUCH.....21
	MANUEL DE L'UTILISATEUR.....41
	MANUAL DEL USUARIO.....61
	MANUAL DO UTILIZADOR.....81
	MANUALE D'USO.....101
	CEBUIKERSHANDLEIDING.....121
	BRUKSANVISNING.....141
	KÄYTTÖOHJE.....161

© 2016 DJO - 157 - Rev A

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




2. MEDICAL BACKGROUND

WHAT IS TENS?

TENS = Transcutaneous Electrical Nerve Stimulation. TENS is an effective method offering Pain Relief. TENS utilises the nervous system's own pain relief mechanisms to ease both acute and chronic pain.

WHEN IS TENS USED?

Good results are observed in acute and long-term pain conditions of many kinds, but above all when the pain originates in joints, the skeleton, muscles, skin, viscera, or nervous system. In some cases, TENS on its own may be sufficient, and in other cases it may provide a valuable complement to other forms of treatment.

HOW DOES IT WORK?

With the Pain Relief programs, impulses are sent through the nerves to block the pain impulses. Pain Relief is usually most efficient during stimulation, but the effect can also last after the program has finished. It also increases the circulation. The Massage programs increase the circulation and decrease the tension in the muscles.

The Pain Relief programs do not always treat the cause of pain. If the pain persists, seek medical advice.

WHEN SHOULD I USE TENS?

The programs in the CEFAR EASY can be used on any occasion when Pain Relief or Muscle Relaxation is needed. All the programs can be used as often as required. Every session should last at least a minimum of 30 minutes up to several hours.

3. GETTING STARTED

PLEASE READ THE PRECAUTIONARY MEASURES - SECTION 12 - BEFORE USING THE CEFAR EASY.

1. INSERT THE BATTERIES

Insert the batteries according to the picture.

2. ATTACH THE ELECTRODES

A. Connect the electrodes to the cable.

B. Attach the electrodes to your body.

C. Connect the cable to the CEFAR EASY.

3. SWITCH THE STIMULATOR ON

Press the ON/OFF button . It can be used for terminating the stimulation at all times, even when the key lock is activated. Always switch off stimulation before removing electrodes from the skin.

4. CHOOSE A PROGRAM

For further information on the programs, see section 5.

5. RAISE THE AMPLITUDE (INTENSITY OF THE STIMULATION)


Push the  button until you reach a comfortable level of stimulation. Always increase the amplitude cautiously. To pause the program, lower the amplitude to 0 with .

BE SURE TO READ THE REST OF THE MANUAL TO FULLY BENEFIT FROM YOUR CEFAR EASY.

4. OVERVIEW OF THE CEFAR EASY

DISPLAY SYMBOLS

 Program selected:
Neck/Shoulder Pain Relief

 Program selected:
Neck/Shoulder Massage

 Program selected:
Lower Back Massage

 Program selected:
Lower Back Pain Relief

 AUTOMATIC KEY LOCK is activated. The key lock is activated 10 seconds after the amplitude has been set. Deactivate by pressing .

 BATTERY CHANGE NEEDED. When you see this symbol in the display, it means that the batteries very soon will run out.

 Broken circuit symbol. See section 9.

 AMPLITUDE. Ongoing amplitude is shown in both digits and as an amplitude bar. The amplitude bar changes to the next square after 5 mA. Lower the amplitude to 0 mA to pause the program.

 Indicates how much is left of the program to run.

5. PROGRAMS

PLEASE READ THE PRECAUTIONARY MEASURES - SECTION 12 - BEFORE USING THE CEFAR EASY.

GB

PAIN RELIEF

LOWER BACK

Program Duration: 30 minutes

Recommended intensity:
Pins-and-needles sensation
without pain. Should not
cause muscle contraction.

**The program can be used
as often as required.**

Electrode Placement and suggested Body Position

PAIN RELIEF

NECK/SHOULDER

Program Duration: 30 minutes

Recommended intensity:
Pins-and-needles sensation
without pain. Should not
cause muscle contraction.

**The program can be used
as often as required.**

Electrode Placement and suggested Body Position

MASSAGE


LOWER BACK

Program Duration:
30 minutes

Recommended intensity:
Muscle vibrations
to light contractions.

**The program can be used
as often as required.**

Electrode Placement and suggested Body Position

MASSAGE

NECK/SHOULDER

Program Duration:
30 minutes

Recommended intensity:
Muscle vibrations
to light contractions.

**The program can be used
as often as required.**

Electrode Placement and suggested Body Position

6. REPLACEMENT OF BATTERIES

When you see the battery symbol  in the display, it means that the batteries very soon will run out. As long as stimulation feels the same as usual you can continue the treatment. It is time to replace the batteries when stimulation starts diminishing or the stimulator stops working.

-
- Switch off the stimulator.
 - Remove the battery cover, take out the batteries.
 - Make sure that the new batteries are placed correctly with respect to the polarity (+ and - symbols). See picture.
 - Dispose the used old batteries at a recycling station.

Rechargeable batteries should be charged when the indicator lights up or the stimulator stops working. NOTE! Never try charging non-rechargeable batteries, due to the risk of explosion.

If the CEFAR EASY is not used for some time (approximately 3 months), the non re-chargeable batteries should be removed from the stimulator.

7. CARE INSTRUCTIONS

Taking care of and cleaning the CEFAR EASY equipment is simple, provided the following instructions are observed:

- Keep stimulator and accessories in the original case when they are not in use. It may, however, be practical to allow the electrodes to remain on the body between treatments. Carbon rubber electrodes can generally remain for 2–3 hours without the electrode gel drying out (does not apply to adhesive gel). They must then be taken off, washed, and dried before being applied again. This is especially important for persons with sensitive skin. In connection with stimulation, make sure that the electrodes are firmly in place.
- When using carbon rubber electrodes, use plenty of electrode gel and avoid drying out by applying tape around all the edges of the electrodes. Rinse the carbon rubber electrodes and the skin with water after use. Do not use detergent for the electrodes.
- Self-adhesive multi-use electrodes are re-moistened if necessary with a few drops of water and kept air-tight (in a plastic bag) on protective paper when they are not in use.
- Never expose the stimulator to water. Wipe it off with a damp cloth if necessary.
- Do not jerk cables or connections.
- The cables are best preserved if left attached to the stimulator between sessions.
- No other maintenance of the device is required. Service life of the device can vary depending on usage conditions. Typical service life is 7 years.
- Never service the device while in use.
- The device should be operated in temperatures between 10°C and 40°C, atmospheric pressures between 50 and 106 kPa, and relative humidity between 30% and 75%.
- The device should be transported and stored in temperatures between -40 °C and 70 °C, atmospheric pressures between 50 and 106 kPa and relative humidity between 10% and 90%.
- Contact manufacturer for assistance in setting up, using or maintaining the equipment or report events

8. ACCESSORIES

Cefar also offers accessories to the CEFAR EASY.
For purchase information, contact your Cefar dealer or visit www.cefar.se.

- **CEFAR EASY BELT**

A lower-back belt is available as an accessory to the CEFAR EASY. The CEFAR EASY belt makes it easier to relieve lower-back pain quickly. The belt is designed to give you warming support and correct placement – every time!

- **CEFAR EASY SELF-ADHESIVE ELECTRODES**

Electrodes get worn out and need to be replaced. It is recommended to replace the electrodes after approximately 20-40 times of usage.

- **Y-CABLE**

With the help of a so-called Y-cable, the number of electrodes can be doubled.

9. TROUBLESHOOTING

THE STIMULATION DOES NOT FEEL THE SAME AS USUAL

- Check that all settings are correct and make sure that the electrodes are correctly placed.
- Slightly change the position of the electrodes.

THE STIMULATION CURRENT FEELS UNPLEASANT

- The skin is irritated.
- The electrodes begin to lose their stickiness and do not stick properly to the skin. Moisten the adhesive surface with a few drops of water before placing on the skin.
- The electrodes are worn out and need to be replaced.
- Slightly change the position of the electrodes.

THE STIMULATION FEELS WEAK OR NOT AT ALL

- Check if the batteries need to be replaced.
- Electrodes are too old and need replacement.

**THE BROKEN CIRCUIT SYMBOL IS SHOWN ON THE DISPLAY SYMBOL**

The broken circuit symbol indicates that the resistance is too high, or that a cable is broken.

- A too high resistance can be caused by a bad connection between the electrodes and your skin, or that the electrodes need to be replaced.
- A cable breakage can be checked by pressing the cable's pins against one another while increasing the amplitude for the corresponding channel to 11 mA. If the amplitude now drops to 0.0 mA and the broken circuit symbol starts flashing, the cable needs to be replaced.

Note! Never increase the amplitude above 20 mA when you check for cable breaks, since this can damage the stimulator.

THE STIMULATOR IS NOT WORKING

If the display to the left appears when you start the stimulator, it means that the stimulator is broken and needs to be replaced. **NOTE!** Do not use the stimulator - contact your Cefar dealer.

10. TECHNICAL DATA

CEFAR EASY is an electrical stimulator intended for Pain Relief and Massage. It delivers a constant current of 60 mA through a load of 1000 Ω .

INFORMATION RELATED TO ELECTROMAGNETIC COMPATIBILITY (EMC).

CEFAR EASY is designed to be used in typical domestic or commercial environments and approved according to the EMC safety standard of EN 60601-1-2.

CEFAR EASY emits very low levels in the radio frequency (RF) interval. Therefore it is not likely to cause any interference in nearby electronic equipment (radios, computers, telephones etc.).

CEFAR EASY is designed to withstand foreseeable disturbances originating from electrostatic discharges, mains supply magnetic fields and radio frequency transmitters (such as mobile telephones).

NUMBER OF CHANNELS	1
CONSTANT CURRENT	Up to a resistance of 1000 Ω . (increased load can reduce the maximum current)
STIMULATION FORMS	Continuous stimulation: modulated pulse duration, modulated frequency.
AMPLITUDE	0-60 mA, in steps of 1 mA (charge maximum: 18 μ C)
WAVEFORM	Symmetrical biphasic square pulse
MAXIMUM FREQUENCY	80 Hz
MAXIMUM PULSE DURATION	300 μ s
POWER SOURCE	2 x 1.5 V AA non-rechargeable or 2 x 1.2 V AA rechargeable batteries
CURRENT CONSUMPTION FOR ONE CHANNEL, 180 μ S, 80 Hz, 30 mA	80 mA
ENVIRONMENT FOR STORAGE, USE AND SHIPPING	Temperature 10 ^o C - 40 ^o C Air humidity 30% - 75% Air pressure 700 hPa-1060 hPa
EXTERNAL DIMENSIONS	52 x 125 x 30 mm
WEIGHT	approximately 130 g incl. batteries
R.M.S. MAX/CHANNEL	10 mA

Warranty: The stimulator is covered by a two-year warranty against manufacturing faults. Subject to normal conditions of use.

11. KEY TO THE STIMULATOR SYMBOLS

Reference number, part number

Follow instruction for uses

Type BF applied parts

Council Directive 2012/19/EU concerning Waste Electrical and Electronic Equipment (WEEE) requires not to dispose of WEEE as municipal waste. Contact your local distributor for information of the unit and accessories

Keep the device dry

Minimum and maximum temperature indications to respect

This equipment complies with all requirements of the Medical Device Directive (93/42/EEC)

Manufacturer name and address, manufacturing year

Power/Pause

12. PRECAUTIONARY MEASURES

12.1. Contraindications

- **Implanted electronic devices.** Do not use the device on patients who have a cardiac pacemaker, implanted defibrillator, or other implanted electronic device, because this may cause electric shock, burns, electrical interference, or death.
- **TENS for Undiagnosed Pain.** Do not use the device as a TENS device on patients whose pain syndromes are undiagnosed.

12.2. Warnings

- **Consult with physician.** Consult with the patient's physician before using the device, because the device may cause lethal rhythm disturbances to the heart in susceptible individuals.
- **Skin condition.** Apply stimulation only to normal, intact, clean, healthy skin.
- **Long term effects.** The long-term effects of chronic electrical stimulation are unknown.
- **Stimulation location.** Stimulation over Neck or Mouth. Do not apply stimulation over the patient's neck (especially the carotid sinus) or the patient's mouth, because this could cause severe muscle spasms resulting in closure of the airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure.
- **Stimulation across Chest.** Do not apply stimulation across the patient's chest, because the introduction of electrical current into the chest may cause rhythm disturbances to the patient's heart, which could be lethal.
- **Across the head.** Since the effects of stimulation of the brain are unknown, stimulation should not be applied across the head, and electrodes should not be placed on opposite sides of the head.
- **Stimulation over compromised skin.** Do not apply stimulation over open wounds or rashes, or over swollen, red, infected, or inflamed areas or skin eruptions (e.g., phlebitis, thrombophlebitis, varicose veins).
- **Stimulation near cancerous lesions.** Do not apply stimulation over, or in proximity to, cancerous lesions.
- **Stimulation over eyes.** Do not apply stimulation directly on the eyes.
- **Environment.** Electronic monitoring equipment. Do not apply stimulation in the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms), which may not operate properly when the electrical stimulation device is in use.
- **Bath or Shower.** Do not apply stimulation when the patient is in the bath or shower. Do not apply stimulation in humid atmosphere exceeding 75% of relative humidity.
- **Sleeping.** Do not apply stimulation while the patient is sleeping.
- **Driving or operating machinery.** Do not apply stimulation while the patient is driving, operating machinery, or during any activity in which electrical stimulation or involuntary muscle contraction can put the patient at risk of injury.

- **Electrosurgical equipment or defibrillators.** Disconnect the stimulation electrodes before using electrosurgical equipment or defibrillators. Otherwise skin burns may be caused below the electrodes and the device might be destroyed.
- **Magnetic Resonance Imaging.** Do not wear electrode or the device during Magnetic Resonance Imaging (MRI) scans as this may result in metal overheating and causing skin burns in the area of the electrode.
- **Flammable or explosive environment.** Do not use the device in areas where there is a risk of fire or explosion, such as oxygen-rich environments, in the vicinity of flammable anaesthetics, etc.
- **Power supply.** Never connect stimulation cables to an external power supply as there is a risk of electric shock.
- **Near other equipment.** Do not use the device beside or stacked on top of any other equipment. If you must use it side by side or on top of another system, you should check that the device works properly in the chosen configuration.
- **Miscellaneous.** Electrodes for Single Patient. Do not share electrodes with other persons. All users should have individual set of electrodes to prevent undesirable skin reactions or disease transmission.
- **Accessories.** Use this device only with the leads, electrodes, and accessories recommended by the manufacturer. Use of other accessories may adversely affect the

performance of the device or may result in stronger electromagnetic emissions or reduce the electromagnetic immunity of the device.

- **No Modification.** No modification of the equipment is allowed.

12.3. Precautions

- **Supervision.** Use this device only under the continued supervision of a licensed practitioner. Electrode placement and stimulation settings should be based on the guidance of the prescribing practitioner.
- **Manufacturer.** The manufacturer does not take any responsibility for any electrode placements other than recommended.
- **Pregnancy.** The safety of electrical stimulation during pregnancy has not been established.
- **Skin irritation.** Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium (gel). The irritation may be reduced by using an alternate conductive medium or alternate electrode placement. Some patients may experience redness under the electrodes after a session. This redness usually disappears within a few hours. Advise the patient to consult the clinician if the skin redness does not disappear after a few hours. Do not start another stimulation session in the same area if the redness is still visible. Don't scratch the redness area.
- **Heart Disease.** Patients with suspected or diagnosed heart disease should follow precautions recommended by their physicians.

- **Epilepsy.** Patients with suspected or diagnosed epilepsy should follow precautions recommended by their physicians.
- **Internal Bleeding.** Use caution when the patient has a tendency to bleed internally, such as following an injury or fracture.
- **After Surgery.** Use caution following recent surgical procedures when stimulation may disrupt the patient's healing process.
- **Over uterus.** Do not place electrodes directly over the uterus or connect pairs of electrodes across the abdomen if you are pregnant. The reason is that, theoretically, the current could affect the foetus's heart (although there are no reports of it being harmful).
- **Lack of sensation.** Use caution if stimulation is applied over areas of skin that lack normal sensation. Don't apply stimulation on patient unable to express themselves.
- **Stimulation.** The stimulator should only be used with skin electrodes intended for nerve and muscle stimulation. Muscle soreness may occur after stimulation but usually disappears within a week.
- **Hot casing or batteries.** Under extreme use conditions, some parts of the casing might reach up to 109 °F (43 °C). Use caution when manipulating the batteries right after device use or when holding the device. There is no particular health risk associated with this temperature besides your comfort.
- **Children.** Keep this device out of the reach of children.
- **Electrode Size.** Do not use electrodes with an active area less than 16 cm², as there will be a risk of suffering a burn injury. Caution should always be exercised with current densities more than 2mA/ cm².
- **Strangulation.** Do not wrap leadwires around your neck, and keep them out of the reach of children. Strangulation may result from entanglement in the leadwires.
- **Tripping.** Care should be used to avoid tripping on lead wires.
- **Damaged Device or Accessories.** Never use the device or any of its accessories if it is damaged (case, cables, etc.) or if the battery compartment is open as there is a risk of electric shock. Carefully inspect the lead wires and connectors prior to each use.
- **Inspect Electrodes.** Inspect electrodes before each use. Replace electrodes when they begin to deteriorate or lose adhesion. Poor contact between the electrodes and the patient's skin increases the risk of skin irritation or burns. Electrodes will last longer if used and stored according to instructions on electrode packaging. Attach the electrodes in such a way that their entire surface is in contact with the skin.
- **Foreign bodies.** Do not allow any foreign bodies (soil, water, metal, etc.) to penetrate the device and the battery compartment.
- **Batteries.** Do not carry batteries in a pocket, purse, or any other place where the terminals could become short-circuited (e.g. by way of paper clip). Intense heat could be generated and injury may result. Never open the battery cover during stimulation

in order to avoid electrical shock. Remove the batteries from the device if you do not intend to use it for a prolonged period of time (more than 3 months). Extended storage of the batteries in the device might lead to batteries and device damage.

- Cable. The cable is best preserved if left attached to the stimulator between sessions. Do not jerk the cable or connection.
- Heat and Cold products. The use of heat or cold producing devices (e.g. electric heating blankets, heating pads or ice packs) may impair performance of the electrode or alter the patient's circulation/sensitivity and increase the risk of injury to the patient.
- Pulled muscles. Do not apply electrodes over pulled muscles. Using the stimulator on a previously extended muscle might further pull such muscle. The higher the stimulation intensity, the higher the risk to further overextends such muscle.
- Additional Precautions for TENS
 - TENS is not effective for pain of central origin, including headache.
 - TENS is not a substitute for pain medications and other pain management therapies.
 - TENS devices have no curative value.
 - TENS is a symptomatic treatment and, as such, suppresses the sensation of pain that would otherwise serve as a protective mechanism.

- Effectiveness of TENS is highly dependent upon patient selection by a practitioner qualified in the management of pain patients.

12.4. Dangers

- Electrodes. Any Electrode with a minimum active area of 16 cm² may be used with this device. Use of an electrode with an area less than 16 cm² can cause burns when the unit is used at higher intensities. Consult your clinician prior to using any electrode less than 16 cm².

12.5. Adverse reactions

- Patients may experience skin irritation and burns beneath the stimulation electrodes applied to the skin.
- Patients may experience headache and other painful sensations during or following the application of electrical stimulation near the eyes and to the head and face.
- Patients should stop using the device and should consult with their physicians if they experience adverse reactions from the device.
- Precaution: Do not disconnect any stimulation cables during a session while the stimulator is switched on. Switch the stimulator off first. Always turn off the stimulator before moving or removing any electrodes during a session.
- Precaution: Do not use electrodes with an active area less than 16 cm², as there will be a risk of suffering a burn injury. Caution should always be exercised with current densities more than 2mA/cm².

- Precaution: Do not apply stimulation in the vicinity of metal. Remove jewelry, body piercings, buckles or any other removable metallic product or device in the area of stimulation. Never use the electrodes contra-laterally, i.e. do not use two pins connected to the same channel on opposite segments of the body.
- Precaution: Never carry out an initial stimulation session on a person who is standing. The first five minutes of stimulation must always be performed on a person who is sitting or lying down. In rare instances, people of a nervous disposition may experience a vasovagal reaction. This reaction is connected with fear of the muscle stimulation as well as surprise at seeing one of their muscles contract without having intentionally contracted it themselves. A vasovagal reaction causes heart to slow down and blood pressure to drop, which can make you feel weak and faint. If this does occur, stop the stimulation and lie down with the legs raised until the feeling of weakness disappears (5 to 10 minutes).
- Precaution: Sudden temperature changes can cause condensation to build up inside the –stimulator. To prevent this, allow it to reach room temperature before use.

14. ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The CEFAR EASY Electrotherapy System is intended for use in the electromagnetic environment defined below. The customer or the user of the CEFAR EASY Electrotherapy System should assure that it is used in such an environment		
Emissions Tests	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Group 1	The CEFAR EASY Electrotherapy System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	CEFAR EASY Electrotherapy System is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes
Harmonic emissions IEC 61000-3-2	Not Applicable - Battery powered	
Voltage fluctuations IEC 61000-3-3	Not Applicable - Battery powered	

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The CEFAR EASY Electrotherapy System is intended for use in the electromagnetic environment specified below. The customer or the user of the CEFAR EASY Electrotherapy System should assure that it is used in such an environment

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Risk assessment on the CEFAR EASY Electrotherapy System indicates the compliance levels claimed are acceptable when ESD-precautionary measures are taken.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	Not Applicable - Battery powered Not Applicable - signal lines less than 3 meters	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+ 1kV differential mode (line to line) + 2kV common mode (line to ground)	Not Applicable - Battery powered	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Not Applicable - Battery powered	Mains power quality should be that of a typical commercial or hospital environment. If the user of the CEFAR EASY Electrotherapy System requires continued operation during power mains interruptions, it is recommended that the CEFAR EASY Electrotherapy System be powered from an uninterrupted power supply.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U_T is the AC mains voltage prior to application of the test level.			

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The CEFAR EASY Electrotherapy System is intended for use in the electromagnetic environment defined below. The customer or the user of the CEFAR EASY Electrotherapy System should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2,5 GHz</p>	<p>[V] V, where $V_i = 3V$</p> <p>[E] V/m, where $E_i = 3V/m$</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the CEFAR EASY Electrotherapy System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = \left[\frac{3.5}{f} \right] \sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = \left[\frac{3.5}{f} \right] \sqrt{P}$ <p>800 MHz to 2,5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CEFAR EASY Electrotherapy System is used exceeds the applicable RF compliance level above, the CEFAR EASY Electrotherapy System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the CEFAR EASY Electrotherapy System.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V] V/m.

Recommended separation distances between portable and mobile RF communications equipment and the CEFAR EASY Electrotherapy System			
<p>The CEFAR EASY Electrotherapy System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the CEFAR EASY Electrotherapy System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the CEFAR EASY Electrotherapy System as recommended below, according to the maximum output power of the communications equipment.</p>			
Rated maximum output power of transmitter P (W)	Separation distance according to frequency of transmitter d (m)		
	150 kHz to 80 MHz $d = \frac{35}{[V_1]} \sqrt{P}$ (where $V_1 = 3V$)	80 MHz to 800 MHz $d = \frac{35}{[E_1]} \sqrt{P}$ (where $E_1 = 3V/m$)	800 MHz to 2,5 GHz $d = \frac{7}{[E_2]} \sqrt{P}$ (where $E_2 = 3V/m$)
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

FCC REQUIREMENTS

Part 15 of the FCC Requirements	
<p>This device complies with Part 15 of the FCC Rules. Operation is subject to the following 2 conditions:</p>	<ul style="list-style-type: none"> • This device may not cause harmful interference • This device must accept any interference received, including the interference that may cause undesired operation.
FCC ID	T9J-RN42
Contains Transmitter Module IC	6514A-RN42

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(9.397 pt)

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Risk assessment on the CEFAR EASY Electrotherapy System indicates the compliance levels claimed are acceptable when ESD-precautionary measures are taken.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	Not Applicable - Battery powered Not Applicable - signal lines less then 3 meters	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+1kV differential mode (line to line) +2kV common mode (line to ground)	Not Applicable - Battery powered	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Not Applicable - Battery powered	Mains power quality should be that of a typical commercial or hospital environment. If the user of the CEFAR EASY Electrotherapy System requires continued operation during power mains interruptions, it is recommended that the CEFAR EASY Electrotherapy System be powered from an uninterrupted power supply.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U_T is the AC mains voltage prior to application of the test level.			

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