

VenaGo™

Pre-Surgery Timeline

for your VenaGo™ portable device to help prevent DVT

- 1** You will receive a call from a DJO® Patient Care Service Representative (PCS rep) ~1 week prior to your surgery. If they leave a voice message, please call them back.
- 2** They will discuss your insurance or financial responsibility. If you have questions, please contact DJO® - 800-452-5065
- 3** Please plug in (and fully charge) the unit prior to surgery to ensure that everything is ready when you return home.

Charging your VenaGo™:

If you wear and operate the VenaGo™ unit while charging expect it to take ~4 hours to fully charge (from a completely depleted battery). If the unit is not in use during charging, expect ~3 hour charge time.

PLEASE NOTE:

Units are not shipped fully charged. Fully charge unit before use.

For Application and Usage Questions

Please Contact Customer Care:

888.405.3251



Copyright © 2021 by DJO, LLC

Individual results may vary. Neither DJO, LLC nor any of its subsidiaries dispense medical advice. The contents of this brochure do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.

MKTIS00-11144 Rev A



T 800.336.6569 F 800.936.6569

1430 Decision Street | Vista, CA 92081-8553 | U.S.A.

DJOGlobal.com/aircast