

AFO LEAF SPRING



AFO Leaf Spring is a pre-fabricated polypropylene ankle-foot orthosis designed to support flaccid drop foot. It provides a semi-rigid section for toe clearance and support, yet has a thin, flexible foot plate for easy trimming. The absence of a heel section makes the Leaf Spring more comfortable to wear and provides a better fit in shoes.

INDICATIONS

Drop foot

FEATURES & BENEFITS

- Injection-molded polypropylene is lightweight
- Variable thickness throughout the orthosis provides strength
- Full length foot plate provides good toe clearance and support
- Thin, flexible foot part is easy to trim with a pair of scissors
- Excellent fit for most types of shoes
- Absence of heel section for improved comfort

ORDERING INFORMATION

| PART # | HEIGHT | LEG | SIZE | FOOT PLATE LENGTH |
|----------|--------|-------|--------------|-------------------|
| O-110127 | 27cm | Right | Small/Medium | 24cm |
| O-110227 | 27cm | Left | Small/Medium | 24cm |
| O-110130 | 30cm | Right | Large | 28cm |
| O-110230 | 30cm | Left | Large | 28cm |