



Polar Care® Glacier™



The Polar Care Glacier Cold Therapy System balances the needs of the patient and the clinician. Its easy to use, feature-rich design makes it a great choice for hospital use. With the addition of ice and water, your patients will enjoy up to 13 hours of cold therapy. Proper use requires an insulation barrier between the pad and the patient's skin. Federal law restricts this device to sale by or on the order of a medical professional.

Common Examples of Use

- Post Operative
- Arthroscopic Procedures
- Reconstructive Procedures
- General Surgery
- Post-Trauma
- Chronic Pain
- Physical Therapy

Features

- Once-a-shift fill requirement increases a hospital's efficiency— 8hr shifts? 12hr shifts? The Glacier has them covered
- Physicians and patients benefit from reduced pain and swelling—follow-up visits focus on recovery, rather than pain management
- Whisper-quiet pump technology allows your patient to rest easy
- Temperature control dial allows for individualized protocol
- Broad selection of premium WrapOn Pads conform to the affected area
- Sterile Polar Dressings available for each part

Category: Devices.

Instructions for Use

Instructions for Use

Polar Care Glacier Patient Instructions


Part Numbers

Part Numbers

Part #	Description
10901	Polar Care Glacier
Combo Units:	PC Glacier with Pad (includes Glacier & 1 single patient use pad)
10902	PC Glacier with Knee/Shoulder Polar Pad Combo
02328	Sterile Polar Dressing, Shoulder
10904	PC Glacier with Multi-Use XL Polar Pad Combo
02348	Sterile Polar Dressing, Multi-Use
10905	PC Glacier with Knee WrapOn Combo
04708	Sterile Polar Dressing, Knee L
10906	PC Glacier with WrapOn Knee Large Pad
04908	Sterile Polar Dressing, Knee L
10907	PC Glacier with WrapOn Knee XL Pad
04908	Sterile Polar Dressing, Knee L
10908	PC Glacier with WrapOn Multi-Use XL Pad
02344	Sterile Polar Dressing, XL Multi-Use
10909	PC Glacier with WrapOn Multi-Use XL Pad with Long Stem
02344	Sterile Polar Dressing, XL Multi-Use
10911	PC Glacier with WrapOn Shoulder Pad
04908	Sterile Polar Dressing, Shoulder L
10912	PC Glacier with WrapOn Shoulder XL Pad
04918	Sterile Polar Dressing, XL Shoulder
10913	PC Glacier with WrapOn Ankle Pad
02508	Sterile Polar Dressing, Ankle
10914	PC Glacier with WrapOn Hip Pad
10680	Sterile Polar Dressing, Knee L
Accessories:	
10698	Replacement Power Supply
10670	Compression strap (60" / 152 cm long)
07640	PC Y Adaptor
09915	4ft Extension Tube
000038-000	PC Bed Hanger

Pad Sizes

WrapOn Polar Pads

Description	Size	Part #
 WrapOn Pad, XL Multi-Use	11.25" x 11.25" / 28 cm x 28 cm	04790

<p><i>Sterile Polar Dressing, XL Multi- Use</i></p>	<p>02344</p>	
<p><i>WrapOn Pad, XL Multi- Use 11.25" x 11.25" / 28 cm x 28 cm Long Stem (16" long)</i></p>	<p>04740</p>	
<p><i>Sterile Polar Dressing, XL Multi- Use</i></p>	<p>02344</p>	
 <p><i>WrapOn Pad, Knee</i></p> <p><i>Sterile Polar Dressing, Knee L</i></p> <p><i>WrapOn Pad, L Knee</i></p> <p><i>Sterile Polar Dressing, Knee L</i></p> <p><i>WrapOn Pad, XL Knee</i></p> <p><i>Sterile Polar Dressing, Shoulder L</i></p>	<p>10.25" x 11.25" / 26 cm x 28 cm</p> <p></p> <p>11.5" x 12" / 29 cm x 30 cm</p> <p></p> <p>12.25" x 19" / 31 cm x 48 cm</p> <p></p>	<p>04700</p> <p>04708</p> <p>04703</p> <p>04708</p> <p>04705</p> <p>04908</p>
 <p><i>WrapOn Pad, Shoulder (Universal)</i></p> <p><i>Sterile Polar Dressing, Shoulder L</i></p> <p><i>WrapOn Pad, Shoulder, w/ Extra Long Straps (9", 12" and 29" / 23 cm, 30 cm and 74 cm long)</i></p> <p><i>Sterile Polar Dressing, Shoulder L</i></p>	<p>10.25" x 11.5" / 26 cm x 29 cm</p> <p></p> <p>10.25" x 11.5" / 26 cm x 29 cm</p> <p></p>	<p>04900</p> <p>04908</p> <p>04903</p> <p>04908</p>

WrapOn

Pad, 13" x 21.5" / 33 cm x 55 cm
XL 10" x 11.75" / 25 cm x 30 cm
Shoulder
(Universal)

04905

Sterile

Polar

Dressing,

XL

Shoulder

04918



WrapOn Pad, Hip

11.5" x 12" / 29 cm x 30 cm

Long strap length = 44" /

112 cm

04750

Short strap length = 26" /

66 cm

Sterile Intelli-Flo Dressing, Hip

10680



WrapOn Pad, Ankle*

8.25" x 9.75" / 22 cm x 25
cm

04730

Sterile Polar Dressing, Ankle

02508



WrapOn Polar Pad, Back**

8.25" x 11" / 21 cm x 28 cm 09800

Sterile Polar Dressing, Back

09810

WrapOn Polar Pad, Back w/ Long Straps (17.5" / 44 cm and 5" / 13 cm Straps)

8.25" x 11" / 21 cm x 28 cm 09805

Sterile Polar Dressing, Back

09810



WrapOn Polar Pad, Hand / Wrist***

8.5" x 9.75" / 22 cm x 25
cm

10750

Sterile
Polar
Dressing,
Small
Ankle

02508

Note: Matching Sterile Polar Dressing listed under each WrapOn Pad. Body part indicated for Sterile Polar Dressing may or may not correlate with body part in pad name, but the sizes have been paired appropriately.

*Not Compatible with Polar Care Cub

**Must be combined with Back Brace PN 1015X or 0733X

***Dressing is included

Polar Pads

	Description	Size	Part #
	Polar Pad, Multi-Use	9.75" x 11.25" / 25 cm x 28 cm	02340
	Sterile Polar Dressing, Multi-Use		02348
	Polar Pad, XL Multi-Use	11.5" x 11.75" / 29 cm x 30 cm	02346
	Sterile Polar Dressing, XL Multi-Use		02344
	Polar Pad, XL Multi-Use, Long Stem	11.5" x 11.75" / 29 cm x 30 cm	02486
	Sterile Polar Dressing, XL Multi-Use		02344
	Polar Pad, Small Rectangle	4.5" x 10.75" / 11 cm x 27 cm	02400

Sterile

Polar

Dressing,

Large

Rectangle

Polar

Pad,

Large

Rectangle

Sterile

Polar

Dressing,

Large

Rectangle

Polar

Pad,

Large

Rectangle

Long

Stem

(16"

long)

Sterile

Polar

Dressing,

Large

Rectangle



Polar Pad, Knee / Shoulder

9.75" x 11.25" / 25 cm x 28 cm

02320

Sterile Polar Dressing, Knee / Shoulder

02328





Polar Pad, Ankle*

8.25" x 9.75" / 21 cm x 25 cm

02500

Sterile Polar Dressing, Ankle

02508

	Polar Pad, Back	6.25" x 11.25" / 16 cm x 28 cm	09900
	<i>Sterile Polar Dressing, Back</i>		09810
	Polar Pad, Eye Mask**	4.25" x 9"	02380
	<i>Sterile Polar Dressing, Eye Mask</i>		02388
	Polar Pad, TMJ**	4.25" x 16" / 16 cm x 28 cm	02460
	<i>Sterile Polar Dressing, TMJ</i>		02468

Note: Matching Sterile Polar Dressing listed under each Polar Pad.

*Not Compatible with Polar Care Cub

**Dressing is included

Resources

Resources

- [Using Cold Therapy to Improve Patient Satisfaction and Reduce Opioid Use Post Surgery](#)