

FORMFIT® PRO ELBOW

Össur Formfit® Pro Elbow provides relief from elbow pain. This 3D compressive elbow sleeve is manufactured using our MotionTech™ technology and features two integrated epicondyle SmartPads™ for more targeted support and cushioning. The RadialStrap™ provides additional compression and support.



MotionTech™



CoolVent™



SmartPad™



FoldLock™

INDICATIONS FOR USE

For elbow conditions that may benefit from compression, such as:

- Medial and lateral epicondylitis
- Elbow pain and swelling


CONTRAINDICATIONS

- Must not be used by individuals for whom compression is contraindicated.

FEATURES AND BENEFITS

- 3D MotionTech™ for compression and dynamic fit keep.
- Epicondyle SmartPads™ for targeted support and cushioning.
- Patented FoldLock™ system for easier and more secure application.
- RadialStrap™ for increased compression and support.

REQUIRED MEASUREMENTS



Measure circumference
12cm (43/4in) below elbow.



Size	Elbow circumference		Blue	Black
	cm	inch	Part#	Part#
XS	17 - 20	6¾ - 7¾	FP301001	FP300001
S	20 - 23	7¾ - 9	FP301002	FP300002
M	23 - 26	9 - 10¼	FP301003	FP300003
L	26 - 29	10¼ - 11½	FP301004	FP300004
XL	29 - 32	11½ - 12½	FP301005	FP300005