

FORMFIT® OA EASE

The Formfit OA Ease is a single upright brace with a Dynamic Force System™ strap designed to relieve pain from mild-to-moderate unicompartmental knee osteoarthritis (OA). Its wraparound design is easy to apply and remove using simple rubber pull-tabs, colour-coded guide dots, and strap tension indicator lines. This low-profile knee brace is comfortable, easy to use, and adapts to a variety of leg anatomies. The Formfit OA Ease can also be machine washed for effortless care.



INDICATIONS FOR USE

Take control of mild-to-moderate OA pain with the Formfit OA Ease, an easy-to-wear knee brace with a simple wraparound design that can be applied and removed with ease.

- Mild-to-moderate unicompartmental knee osteoarthritis and/or degenerative meniscal tears
- Other conditions that may benefit from unicompartmental unloading

FEATURES AND BENEFITS

- Easy Wraparound Design - Simple rubber pull-tabs, colour-coded guide dots and a convenient hinge alignment guide make application and removal simple
- Relieves OA Pain - 3-Point Leverage System with Dynamic Force System™ strap and shells unload the affected knee compartment
- Comfortable - Lightweight and breathable fabric feels great against the skin and can be machine washed for effortless care
- Discreet – Low-profile brace with a slim hinge can be worn under clothes
- Control – Fine-tune strap tension for individual pain relief, with marked dosing indicator lines
- Extension kits and liners are available to further enhance fit and comfort

REQUIRED MEASUREMENTS



Measure calf circumference 15cm below mid-patella. Leg in extension.



Measure thigh circumference 15cm above mid-patella. Leg in extension.

FORMFIT® OA EASE

ORDERING INFORMATION

Size	Calf circ. 15cm below mid patella	Thigh circ. 15cm above mid patella	Right Medial	Left Medial	Right Lateral	Left Lateral
XS	28 - 32	36 - 41	FOA0518711	FOA0618711	FOA0528711	FOA0628711
S*	32 - 36	41 - 46	FOA0518712	FOA0618712	FOA0528712	FOA0628712
M	36 - 41	46 - 53	FOA0518713	FOA0618713	FOA0528713	FOA0628713
L*	41 - 46	53 - 60	FOA0518714	FOA0618714	FOA0528714	FOA0628714
XL	46 - 51	60 - 67	FOA0518715	FOA0618715	FOA0528715	FOA0628715
XXL*	51 - 56	67 - 75	FOA0518716	FOA0618716	FOA0528716	FOA0628716

*All sizes are available but for consolidated stock option this size can be achieved by combining the smaller size with accessories below.

ACCESSORIES

Description	Part#
ActiveGrip Calf Liner Kit 2MM	FOA000101
ActiveGrip Calf Liner Kit 6MM	FOA000102
Doeskin Thigh Liner Kit 2MM	FOA000201
Doeskin Thigh Liner Kit 6MM	FOA000202
Small Extension Kit 2 Pieces	FOA000301
Small Extension Kit 10 Pieces	FOA000302
Large Extension Kit XS-SM (Right Medial/Left Lateral)	FOA000401
Large Extension Kit MD-LG (Right Medial/Left Lateral)	FOA000402
Large Extension Kit XL-XXXL (Right Medial/Left Lateral)	FOA000403
Large Extension Kit XS-SM (Right Lateral/Left Medial)	FOA000404
Large Extension Kit MD-LG (Right Lateral/Left Medial)	FOA000405
Large Extension Kit XL-XXXL (Right Lateral/Left Medial)	FOA000406
Wash Bag	WB000001
Suspension Strap	OAW00007
Flexion/Extension Stop Kit	B-301040047

REPLACEMENT PARTS

Description	Part#
Alligator Clip (Right Medial/Left Lateral)	FOW00501
Alligator Clip (Right Lateral/Left Medial)	FOW00502
Hinge Pad	B-304120101