


LumboLoc® Support

(LumboLoc® Forte ES)

Increased stabilization for optimal spinal posture and relief of the lumbosacral spine



Integrated dorsal corset stays straighten the lumbar spine




Lumbosacral Pad



Abdominal pad



Functional strap



Stabilization and relief

INDICATIONS

- General pain in the lumbar spine region (acute and chronic)
- Moderate lumbar sciatica with slight muscular deficiencies in cases of intervertebral disk protrusions / herniated disks
- Moderate (pseudo)radicular lumbar syndrome
- Moderate degeneration / muscular insufficiency
- Condition following herniation, e. g. post-operative
- Moderate facet syndrome
- Moderate lumbar deformity
- Ligamento-tendomyopathy
- Spondylolysis / spondylolisthesis, grade 1

BENEFITS AND MODES OF ACTION

- Integrated corset stays straighten the lumbar spine
- Increased stabilization thanks to individually adjustable tensioning straps
- Additional stabilization and support thanks to the lumbosacral pad
- Greater straightening thanks to increased abdominal pressure exerted by the abdominal pad
- Height-adjustable functional strap for indication-specific pressure distribution
- Practical finger loops and strap system allows the orthosis to be fitted easily

ACCESSORIES

Extension panel: 12984881081000

- If the biggest size available is too small, the **extension panel** allows to extend the maximum torso circumference of the back brace by around 10 centimeters.

SIZE	1	2	3	4	5	6	ORDER NO.
Circumference trunk (cm)	70 – 80	80 – 90	90 – 100	100 – 110	110 – 120	120 – 130	121 8399 0080 00°

° = Size (1 – 6)

