

FORMFIT® PRO KNEE OA

Össur Formfit® Pro Knee OA is a unique, compressive knee sleeve with a Dynamic Force Strap system, designed to provide pain relief from the early stages of osteoarthritis (OA). The dynamic, 3-D MotionTech™ knit is supportive and meets certified medical grade compression requirements, and has been proven to be superior to the leading brands at moisture wicking and breathability, keeping your knee cool during wear. CoolVent™ technology in the popliteal area enhances breathability and comfort, while the distinct patella graphics help with intuitive placing of the sleeve. The integrated StayBilizer™ hinge, constructed from lightweight, flexible EVA material, is suitable for sports, such as football.



MotionTech™



CoolVent™



StayBilizer™



FoldLock™

INDICATIONS FOR USE

For osteoarthritis and other knee conditions which may benefit from unloading and/or compression, such as:

- Medial or lateral knee pain
- Early or mild unicompartmental osteoarthritis of the knee
- Symptomatic, unicompartmental meniscus degeneration

CONTRAINDICATIONS

- Must not be used by individuals for whom compression is contraindicated.

FEATURES AND BENEFITS

- Patented DFS system.
- Certified medical grade compression (20-30 mmHg).
- Lightweight CoolVent™ knit in popliteal area for added breathability and comfort.
- Patented FoldLock™ system for easier and more secure application.
- Sports-compatible StayBilizer™ hinge made from light, flexible EVA material for optimum shape and fit.
- Patella graphics for intuitive placement.
- Available in 5 standard sizes and 3 comfort sizes.

REQUIRED MEASUREMENTS



Measure calf circumference 15 cm (6") below mid-patella. Leg in extension.



Measure thigh circumference 15 cm (6") above mid-patella. Leg in extension.

Size	Calf circumference 15 cm (6") below mid patella		Thigh circumference 15 cm (6") above mid patella		Right Medial, Left Lateral	Right Lateral, Left Medial
	cm	inch	cm	inch		
XS	26 - 31	10¼ - 12¼	36 - 41	14□ - 16□	FP110401	FP110301
S	31 - 36	12¼ - 14□	41 - 46	16□ - 18□	FP110402	FP110302
M	36 - 41	14□ - 16□	46 - 51	18□ - 20□	FP110403	FP110303
M Comfort	36 - 41	14□ - 16□	51 - 56	20□ - 22	FP110413	FP110313
L	41 - 46	16□ - 18□	51 - 56	20□ - 22	FP110404	FP110304
L Comfort	41 - 46	16□ - 18□	56 - 60	22 - 23□	FP110414	FP110314
XL	46 - 50	18□ - 19¾	56 - 60	22 - 23□	FP110405	FP110305
XL Comfort	46 - 50	18□ - 19¾	60 - 64	23□ - 25¼	FP110415	FP110315