




Quad Canes

ProBasics Quad Canes feature a unique four-point design which allows for greater weight-bearing and stability than a single-point cane. The handle rotates allowing for a proper fit with right or left hand use. Small and large base models available in black or silver finish with a limited lifetime warranty.



Easy one-button height adjustment



Contoured vinyl hand grip

FEATURES

- Unique four-point design allows for greater weight-bearing and stability than a single-point cane
- Non-skid rubber tips on all four feet enhance durability for long term use
- Small base design allows for easier maneuverability while the large base design provides more stability
- Handle rotates for left or right hand use
- Quad canes stand up on their own when not in use
- Easy push button height adjustment in 1-inch increments with a locking, anti-rattle ring allowing for quick, safe adjustment
- Comfortable, contoured vinyl handle with sculpted grip reduces chance of slippage
- Lightweight aluminum upper shaft with welded steel base
- Tool-free left to right hand adjustments
- Bariatric model features a double-plated base for increased stability and support for patients up to 500 pounds
- HCPCS* Code: E0105

SPECIFICATIONS

Handle Height Adjustment Range:	Small & Large Base: 29" – 38" Bariatric: 27" – 36"
Cane Tip:	Small & Large Base: 3/4" Bariatric: 3/4"
Weight Capacity:	Adult: 300 lb Bariatric: 500 lb
Product Weight:	Small Base: 2.3 lb Large Base: 2.7 lb Bariatric Small Base: 3.2 lb Bariatric Large Base: 3.6 lb
Carton Dimensions:	Small Base: 25.6" x 10.6" x 8.7" Large Base: 25.6" x 13.8" x 11.8" Bariatric Small Base: 24.4" x 16.5" x 6.3" Bariatric Large Base: 24.2" x 12.2" x 8.9"
Shipping Weight:	Small Base: 14.6 lb Large Base: 17.8 lb Bariatric Small Base: 14 lb Bariatric Large Base: 8.5 lb
Warranty:	Limited Lifetime on frame 90 days on all parts

Item #	Description	Color	UPC	UOM
CNQSB	ProBasics Small Base Quad Cane	Black	815067071180	6/cs
CNQLB	ProBasics Large Base Quad Cane	Black	815067071197	6/cs
CNQSBS	ProBasics Small Base Quad Cane	Silver	815067071203	6/cs
CNQLBS	ProBasics Large Base Quad Cane	Silver	815067071210	6/cs
CNQHSBB	ProBasics Bariatric Small Base Quad Cane	Black	815067071227	4/cs
CNQHLBB	ProBasics Bariatric Large Base Quad Cane	Black	815067071234	2/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.