



Genesis Full Shell Walker



Sizes: xs, s, m, l, xl

The Genesis Walker is the lightest full shell pneumatic walker on the market. It has been computer engineered to provide strength where required while eliminating extra weight. Similar to popular athletic shoes, the Genesis Walker uses a dynamic air chamber in the heel of the tread which actively absorbs the shock of heel strike and allows for a more comfortable walking experience.

Also Available in a Mid-Calf Version

Common Examples of Use

- Ankle sprains and strains
- Forefoot sprains
- Stable foot fractures
- Distal tibial fractures

Product Features

- Built-in air pump
- Durable hard shell design
- Trim-able shell
- Narrow heel width
- Dynamic air stride technology
- Dual air chambers
- Scalloped no-slip tread

Category: Walker Boots.

Instructions for Use

[Instructions for Use](#)

Part Numbers

Part Numbers

Part #	Description
BL5150XX	Genesis 3-Strap, XS-XL
BL5140XX	Genesis 4-Strap, XS-XL
100507-010	Genesis Cool Liner - Tall, XS
100507-020	Genesis Cool Liner - Tall, S
100507-030	Genesis Cool Liner - Tall, M
100507-040	Genesis Cool Liner - Tall, L
100507-050	Genesis Cool Liner - Tall, XL

Literature

[Genesis Walker Brochure](#)

Product Sizing Chart

Genesis Full Shell Walker Sizing Chart			
Size #	Size	Men	Women
XX = 01	XS	1 – 2.5	1 – 3.5
XX = 03	S	2.5 – 5.5	3.5 – 6
XX = 05	M	5.5 – 9	6 – 9.5
XX = 07	L	9 – 12.5	9.5 – 13
XX = 09	XL	12.5 – 17	13 +

Measurements according to shoe size