



Hi-Performance Knit Knee Support



Sizes: xs, s, m, l, xl, xxl

The Hi-Performance Knit Support provides compression, warmth, and patella stabilization. It features three-dimensional knitting for a breathable, comfortable fit and a contoured silicone patella donut for support.

Common Examples of Use

- Mild knee sprains
- Chondromalacia
- Mild patellofemoral osteoarthritis
- Osgood Schlatters disease

Product Features

- Breathable 3-dimensional knit
- Silicone patella donut
- Medial and lateral stays to help prevent migration

Category: Soft Supports.

Instructions for Use

[Instructions for Use](#)

Part Numbers

Part Numbers

Part #	Description
2804X	Hi-Performance Knit Support

Product Sizing Chart

Hi-Performance Knit Support			
Size #	Size	Thigh Circumference	
X = 1	XS	15" – 16 1/8"	(38 – 41 cm)
X = 2	S	16 1/8" – 17 3/8"	(41 – 44 cm)
X = 3	M	17 3/8" – 18 1/2"	(44 – 47 cm)
X = 4	L	18 1/2" – 19 3/4"	(47 – 50 cm)
X = 5	XL	19 3/4" – 20 7/8"	(50 – 53 cm)
X = 6	XXL	20 7/8" – 22"	(53 – 56 cm)
<i>Circumference taken 5.5" (14 cm) above mid patella.</i>			