



Roadrunner™ Soft Knee Brace



Sizes: xs, s, m, m+, l, xl, xxl

The Roadrunner provides control and support for patients with ligament injuries and instabilities and is frequently used for post-op treatment and mild OA. Measures 16" in length.

Common Examples of Use

- Most ACL, PCL, MCL, LCL injuries and instabilities
- Meniscal injuries
- Mild OA

Product Features

- Available in Airmesh® and neoprene
- Polycentric hinges with adjustable flexion extension stops
- Sleeve and wraparound versions available
- Standard and Post-Op hinge options

Category: Adjustable ROM.

Instructions for Use

[Instructions for Use](#)

Part Numbers

Part Numbers

Part #	Description
1411X	Roadrunner, Airmesh Pull-on
1414X	Roadrunner, Airmesh, Open Back Pull-on
1415X	Roadrunner, Airmesh, Open Back, Wraparound
0694X	Roadrunner, Neoprene Pull-on
0696X	Roadrunner, Neoprene, Open Back, Pull-on
0697X	Roadrunner, Neoprene, Open Back w/ Adjustable Horseshoe Pull-on
0700X	Roadrunner, Neoprene, Open Back, Wraparound
2182X	Roadrunner, Neoprene, Patella Stabilizer, Open Back, Pull-on
73562	Hinge Cover, Airmesh
73561	Hinge Cover, Neoprene

Product Sizing Chart

Roadrunner Size Chart		
Size #	Size	Thigh Circumference
X = 1	XS	12" - 15" (30 - 38 cm)
X = 2	S	15" - 18" (38 - 46 cm)
X = 3	M	18" - 21" (46 - 53 cm)
X = 4	L	21" - 24" (53 - 61 cm)
X = 5	XL	24" - 27" (61 - 69 cm)
X = 6	XXL	27" - 30" (69 - 79 cm)

Circumference taken 6" (15 cm) above mid patella.

Brace measures 16" in length.