

# AYRE

## WIRELESS AIR COMPRESSION RECOVERY BOOTS



Feel relief aid anywhere, anytime with **Compex Ayre** compression boots. Designed for elite and recreational athletes on-the-go, trainers and coaches. **Compex Ayre** uses a dynamic pulsing compression up the limbs in waves for the duration of its use, which helps promote quick healing.

USE TO HELP TREAT / PREVENT: Sore or aching muscles, pain and swelling of the muscles.

BEST FOR: Aiding in the relief of sore or aching muscles from exercise or recreational activities helping you to be ready to go at full speed the next day.



### FEATURES

- Ultra portable, no tethers to a power outlet during compression
- 4 Separate compression chambers can be inflated separately or simultaneously
- 2 modes: Inflate the four chambers one by one or all four chambers at the same time
- Lithium ion rechargeable battery for up to 60 minutes of recovery
- High quality exterior and interior linings make cleaning easy, and help keep unwanted odors in check
- Integrated compression controls remove the bulk of a separate, wired unit
- Powerful compression: The pressure can be adjusted by increments of 10 from 0mmHg up to a max of 120mmHg

### 2 YEARS WARRANTY

The device may be used while it is charging

BOX CONTAINS: COMPRESSION BOOTS [Quantity: 2]  
CHARGER · CASE · QUICK START GUIDE

SIZE: 44.5x39.5x37 cms

WEIGHT: 2.2 Kgs.

PN	DESCRIPTION	EAN13 (GTIN-13)	GS1-128	GTIN-14
INTL-CX201FB01	COMPEX COMPRESSION BOOTS SET INTL VERSION	0190446701439	(01) 00190446701439	00190446701439

Not medical device