



MARVEL

Before DonJoy® became a global leader in sports medicine, there was Mark Nordquist, a line captain with the Philadelphia Eagles, and Ken Reed, a local lawyer, who started assembling the neoprene sleeves that would serve as the company's first knee, ankle, and elbow braces. Used by elite athletes worldwide to enhance performance, protect the body, and aid recovery, DonJoy's comprehensive range of products can be seen anywhere from the NCAA to the NFL to the global stage of the Olympics. And now, DonJoy® Advantage brings these leading sports medicine technologies direct to the consumer.

For the first time, DonJoy Advantage is expanding its line of bracing and supports to youth through the Marvel franchise. Every year, millions of children are affected with sport-related injuries, such as sprains and strains, ligament tears, and broken bones. The need for preventative, protective, and recovery products is addressed by designing seven unique products that feature Marvel's signature characters.



PATELLA KNEE SLEEVE

Neoprene knee sleeve with patella donut for mild joint support



RECOMMENDED FOR: Strains, Tendonitis, Patella Support


FEATURES & BENEFITS

- 3mm Neoprene aids in a compressive fit
- Pull tabs for easy on/off
- Mesh popliteal for a comfortable fit
- Calf expansion panel aids in flexibility

SUPPORT:

MILD
MODERATE
MAXIMUM

PART NUMBER	SIZE	COLORWAY
DA191KS01-SPI-PED	PEDIATRIC	SPIDERMAN
DA191KS01-CAP-PED	PEDIATRIC	CAPTAIN AMERICA
DA191KS01-SPI-Y	YOUTH	SPIDERMAN
DA191KS01-CAP-Y	YOUTH	CAPTAIN AMERICA

	PEDIATRIC	YOUTH
	11 - 12"	12 - 13"
	28 - 30.5 cm	30.5 - 33 cm

*Size based on knee circumference

ELASTIC KNEE

Lightweight knit sleeve designed to provide compression and mild support.



RECOMMENDED FOR: Sprains and Strains, Swelling, Soreness, Arthritis

FEATURES & BENEFITS

- Stretch elastic fabric designed for lightweight and breathable support
- Expansion panels aid in freedom of movement
- Universal design, fits the left or right knee


SUPPORT:

MILD

MODERATE

MAXIMUM

PART NUMBER	SIZE	COLORWAY
DA161KS01-IRM-PED	PEDIATRIC	IRON MAN
DA161KS01-HLK-PED	PEDIATRIC	HULK
DA161KS01-IRM-Y	YOUTH	IRON MAN
DA161KS01-HLK-Y	YOUTH	HULK

	PEDIATRIC	YOUTH
	11 - 12"	12 - 13"
	28 - 30.5 cm	30.5 - 33 cm

*Size based on knee circumference

COMFORT WRIST BRACE

Wrist brace with cushion foam lining and removable aluminum palmar stay helps support the wrist.




RECOMMENDED FOR: Carpal Tunnel, Sprains and Strains, Tendonitis, Wrist Pain

SUPPORT: MILD MODERATE MAXIMUM

FEATURES & BENEFITS

- Removable aluminum palmar stay aids in immobilization
- Cushion foam lining helps provide a comfortable fit
- 3 trimmable straps allow for adjustability
- Knit expansion panel provides added flexibility and comfort
- Easy to apply, slip-on design
- Anatomical design, right and left specific

RIGHT	LEFT	SIZE	COLORWAY
DA191WB01-CAP-PED-R	DA191WB01-CAP-PED-L	PEDIATRIC	CAPTAIN AMERICA
DA191WB01-CAP-PED-R	DA191WB01-SPI-PED-L	PEDIATRIC	SPIDERMAN
DA191WB01-CAP-Y-R	DA191WB01-CAP-Y-L	YOUTH	CAPTAIN AMERICA
DA191WB01-SPI-Y-R	DA191WB01-SPI-Y-L	YOUTH	SPIDERMAN

	PEDIATRIC	YOUTH
	3.5 - 4.5"	4.5 - 5.5"
	9 - 11 cm	11 - 14 cm

*Size based on wrist circumference

FIGURE-8 ANKLE SUPPORT

Compressive sleeve and figure-8 strap design helps stabilize and support the ankle.



RECOMMENDED FOR: Sprains and Strains, Swelling, Tendonitis

FEATURES & BENEFITS

- Mesh compression sleeve for breathability and comfort
- Figure-8 straps aid in ankle support
- Open heel design for anatomical fit
- Universal design, fits the left or right ankle

SUPPORT:

MILD

MODERATE

MAXIMUM

PART NUMBER	SIZE	COLORWAY
DA161AV03-SPI-PED	PEDIATRIC	SPIDERMAN
DA161AV03-CAP-PED	PEDIATRIC	CAPTAIN AMERICA
DA161AV03-SPI-Y	YOUTH	SPIDERMAN
DA161AV03-CAP-Y	YOUTH	CAPTAIN AMERICA

	PEDIATRIC	YOUTH
	5.75 - 6.75"	6.75 - 7.75"
	15 - 17 cm	17 - 19 cm

*Size based on ankle circumference

YOUTH ARM SLING

Arm sling with lightweight envelope and strap release clip designed for easy application and comfortable fit.



RECOMMENDED FOR: Immobilization and support of the arm, shoulder and elbow joints

FEATURES & BENEFITS

- Strap release clip aids in on/off application
- Nylon, Spandex envelope helps with a comfortable fit
- Adjustable strap length aids in a proper fit
- Universal design fits right or left arm

SUPPORT: MILD MODERATE MAXIMUM

PART NUMBER	SIZE	COLORWAY
DA191AR01-AVE-PED	PEDIATRIC	AVENGERS
DA191AR01-HLK-PED	PEDIATRIC	HULK
DA191AR01-CAP-PED	PEDIATRIC	CAPTAIN AMERICA
DA191AR01-IRM-PED	PEDIATRIC	IRON MAN
DA191AR01-AVE-Y	YOUTH	AVENGERS
DA191AR01-HLK-Y	YOUTH	HULK
DA191AR01-CAP-Y	YOUTH	CAPTAIN AMERICA
DA191AR01-IRM-Y	YOUTH	IRON MAN

	PEDIATRIC	YOUTH
	Up to 9"	Up to 11"
	Up to 23 cm	Up to 28 cm

*Size based on arm length

REUSABLE COLD PACK (ROUND)

Reusable cold pack remains flexible when frozen and helps to provide soothing relief to pain and swelling associated with injury.



RECOMMENDED FOR: Muscle, Joint Pain and Swelling

FEATURES & BENEFITS

- Cold therapy designed to provide relief to pain and swelling
- Malleable cold pack when frozen designed to conform around anatomy during use
- Reusable cold pack designed to be refrozen after each use

PART NUMBER	COLORWAY
DA193CT01-SPI	SPIDERMAN
DA193CT01-CAP	CAPTAIN AMERICA
DA193CT01-IRM	IRON MAN
DA193CT01-HLK	HULK

REUSABLE COLD PACK (SMALL)

Reusable cold pack remains flexible when frozen and helps to provide soothing relief to pain and swelling associated with injury.



RECOMMENDED FOR: Muscle, Joint Pain and Swelling

FEATURES & BENEFITS

- Cold therapy designed to provide relief to pain and swelling
- Malleable cold pack when frozen designed to conform around anatomy during use
- Reusable cold pack designed to be refrozen after each use

PART NUMBER
DA193CT02-AVE

KNEE PRODUCTS



ANKLE / WRIST PRODUCTS



ARM PRODUCTS



COLD PACK PRODUCTS



PACKAGING



DJOGLOBAL.COM | MKT00-10733 REV A

Individual results may vary. Neither DJO, LLC, nor any of its subsidiaries dispense medical advice. The contents of this catalog do not constitute medical, legal, or any other professional advice. Information related to various health, medical, and fitness conditions and their treatment is not meant to be a substitute for the advice provided by a physician or other medical professional. You should not use the information contained herein for diagnosing a health or fitness problem or disease. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.

CONTACT:
 DJO Canada
 6485 Kennedy Road,
 Mississauga, ON, L5T 2W4
 1-866-866-5031



©MARVEL