

NEW  
DESIGN

# AKS™ Knee Support with Metal Hinges & Straps



The AKS Knee Support offers your patients the unique ability to adjust compression and fit from top to bottom.



## Adjustable Knee Sleeve

The distal and proximal openings allow for easy application and ability to fit various patient shapes. In addition, compression can be adjusted to the desired level.

## Adjustable Patellar Buttress

Soft horseshoe buttress is repositionable and works in concert with central sleeve region to provide circumferential compression around the knee and focal compression around the patella.

## Medial and Lateral Metal Hinges & Straps

Metal hinges and outer closure straps provide medial and lateral stability. Internal Skinloc™ pad and metal hinges function synergistically to resist distal migration and vertical bunching. Condylar pads enhance comfort.



## Indications:

- Arthritis of the knee
- Bursitis of the knee
- Chondromalacia patella
- Lateral or medial collateral ligament sprains
- Lateral or medial instability of the knee
- Patellar tendinitis

Bilateral design fits left or right knee.

Neoprene version provides greater compression and warmth.

CoolFlex version is breathable and lightweight for superior comfort.

Neoprene	Size	Knee Circ.*	CoolFlex
117111	x-small	10" - 12"	117121
117112	small	12" - 14"	117122
117114	medium	14" - 16"	117124
117115	large	16" - 18"	117125
117116	x-large	18" - 20"	117126
117117	xx-large	20" - 22"	117127
117118	xxx-large	over 22"	117128

Suggested HCPCS: L1820

CE  
rev 12/10



\*measure circumference of knee at mid patella

Medical Specialties, Inc.  
4600-K Lebanon Rd.  
Charlotte, NC 28227  
p: 800-334-4143 f: 800-694-9060  
p: 704-694-2434 f: 704-694-9060  
e-mail: request@medspec.com  
www.medspec.com  
Patent pending

©2010 Medical Specialties, Inc.