



# Gripper™ 12" Hinged Knee Brace



The Gripper™ Hinged Knee Brace provides your athletes and patients with the medial and lateral stability they need while resisting distal migration.



### Compliance & Migration

Internal SkinLoc™ pad works in concert with the outer closure strap to grip the proximal end of the calf muscle and inhibit the brace from migrating distally, thereby enhancing patient compliance

### Durability

Straps are made of high tension elastic that allows for numerous applications without fraying or losing closure strength. Metal hinge pockets are reinforced on the inside with ballistic nylon material at wear points to prevent hinge from wearing through the pocket.

### Ease of Application & Comfort

The wraparound anterior closure allows for ease of application and an optimum fit. Open popliteal minimizes material bunching during flexion to improve comfort.



The Gripper 12" is available in two different base materials: Neoprene and CoolFlex™. The neoprene version provides warmth and superior compression. The Gripper 2 is made with a thinner and lighter neoprene. The CoolFlex version is breathable and hypoallergenic.

### Indications:

- LCL or MCL sprains
- Lateral or medial instability
- Limit hyperextension

Suggested HCPCS: L1820

Gripper Neoprene 3/16"	Gripper CoolFlex	Size	Knee Circ.	Gripper 2 Neoprene 1/8"
117011	117031	X-Small	10" - 12"	117151
117012	117032	Small	12" - 14"	117152
117014	117034	Medium	14" - 16"	117154
117015	117035	Large	16" - 18"	117155
117016	117036	X-Large	18" - 20"	117156
117017	117037	XX-Large	20" - 22"	117157
117018	117038	XXX-Large	22" - 24"	117158
117019	117039	4X-Large	24" - 26"	117159
117010	117030	5X-Large	26" - 29"	

\*measure circumference of knee at mid patella

