



# Wrist Lacer™ II

## wrist support



With its patent pending design, the new Wrist Lacer II provides equal tension across all laces offering supreme comfort, immobilization and ease of application.

### Adjustability & Immobilization

Blue palmar pad offers enhanced support by matching the contour of the palm. The dorsal and palmar aluminum stays are removable to allow bending by medical professionals to create a custom fit.

### Lacer II Closure System

Unique lace closure design provides equal tension across all laces to easily generate uniform compression for unparalleled support. The tension is easily adjusted for patients with unusual forearm and wrist circumferences.

### Optimum Comfort

Pad inside stockinette protects tendons by cushioning distal lace at base of thumb. The polypropylene felt liner and perforated suede allow moisture to evaporate away from the skin thereby enhancing patient compliance.



8" wrist support



10 1/2" wrist & forearm support

### Indications:

- Carpal tunnel syndrome
- Wrist sprains or strains
- Post cast removal
- Postoperative care

Suggested HCPCS: L3908

Size	8" Back Left	8" Black Right
Pediatric	223329	223339
X-Small	223321	223331
Small	223322	223332
Medium	223324	223334
Large	223325	223335
X-Large	223326	223336
XX-Large	223327	223337

Sizing Chart	
Size	Wrist Circ.
xx-small	under 4 3/4"
x-small	4 3/4" - 5 1/2"
small	5 1/2" - 6 1/4"
medium	6 1/4" - 7"
large	7" - 7 3/4"
x-large	7 3/4" - 8 1/2"

Size	10 1/2" Black Left	10 1/2" Black Right
small	223362	223372
medium	223364	223374
large	223365	223375
x-large	223366	223376



Medical Specialties, Inc. 4600-K Lebanon Rd. Charlotte, NC 28227 p: 800-334-4143 f: 800-694-9060 p: 704-694-2434 f: 704-694-9060

e-mail: request@medspec.com www.medspec.com Patent Pending

© 2015-2017 Medical Specialties, Inc.

