

# 2019-2020 PRODUCT GUIDE

Hand		RM Restorative Hand	RM Prosperity Hand	NeuroFlex® Restorative Flex Hand	NeuroFlex® Restorative ThumbEase	RM Cock-Up Wrist Splint	RM Cock-Up Hand Splint with Finger Extension Roll	RM BendEase Hand	RM Pediatric HyperHand	Restorative Wrist
Item #		#20736	#20732	#20721 (adult) #20720 (pediatric)	#20737 / #20740	#20739	#20739roll	#20715	#20717	RM7001
Wrist/ Hand/ Fingers	Mild/ Moderate Contracture	✓	✓	✓	✓	✓	✓	✓	✓	
	Severe Contracture	Remold. Thumb determines splint choice	Remold	Moderate to Closed Palm	Remold		Remold	Remold by Hand	Remold	Readjust Bolt
Pediatric Sizes Available				✓					✓	
Finger	Severe Flexion in 4th & 5th Digit		✓	✓						
	Severe Thumb ADD			✓	✓	✓	✓		Remold	
Moldable By Hand								✓		
Motor or Movement Disorders Related To:										
Cerebral Palsy		✓	✓	✓	✓	✓	✓	✓	✓	
Multiple Sclerosis		✓	✓	✓	✓	✓	✓	✓	✓	
Arthritis		✓	✓	✓	✓	✓	✓	✓	✓	✓
CVA or Stroke		✓	✓	✓	✓	✓	✓	✓	✓	
Carpel Tunnel Syndrome		Remold Wrist	Remold Wrist	Remold Wrist	Remold Wrist	Molded to 30 degrees ext.	Molded to 30 degrees ext.	Molded to 30 degrees ext.	Remold Wrist	

\* This chart is intended to provide caregivers a snapshot of different product features. It may not be comprehensive and should not be used as the sole factor in product selection and use. Each individual's needs are unique.