

Restorative™ Hand

Effectively and comfortably treat any wrist/hand/finger lost range and tension due to traumatic brain injury, stroke, spinal cord injury, cerebral palsy, multiple sclerosis, or Parkinson's disease.

Features:

- Flex design encourages joints to relax.
- Heat mold to provide mild stretch on the joints, then remold to gradually correct joints to normal alignment.
- Provides prolonged low load passive stretch to relengthen shortened tissue.
- Kydex base extends to end of thumb.
- Comfortably supports and treats osteo and rheumatoid arthritic joints.
- Prevents deformities if treated in early stages.
- Includes finger seperators (on adult models).
- Colour coded straps ease in donning.



Adult Model



Pediatric Model

Now available with dorsal attachment in all sizes, ideal for Constrained Induced Movement Therapy.



Dorsal Attachment


Product Codes:

Adult: 20736- (size)(left/right)

Pediatric: 20717 - (size)(left/right)

Adult w/ dorsal 20736- (size)(left/right) - MOD

Pediatric w/dorsal: 20717 - (size)(left/right) - MOD



Adult Sizing

Size	Width of MCP	Wrist to Fingers
Small	Up to 3"	Up to 6.5"
Medium	Up to 3.5"	Up to 7"
Large	Up to 3.75"	Up to 8"

Both Adult and Pediatric Sizes are Available in Left or Right

Pediatric Sizing

Size	Width of MCP	Wrist to Finger
Infant	1.75" to 2"	Up to 2.5"
Toddler Jr.	1.75" to 2"	Up to 3"
Toddler	2" to 2.25"	Up to 3.5"
Child	2.25" to 2.5"	Up to 4"
Youth	2.75" to 3"	Up to 5"

Phone: 1-866-334-4441
Fax: 1-888-484-4431



Online: www.poportho.com
Email: csr@poportho.com

Professional Orthopedic Products