

## REBOUND® DUAL ST

The Rebound® DUAL ST offers simplicity without compromise. With its malleable yet supportive frame, Rebound DUAL ST is the go-to for knee ligament support.



### INDICATIONS FOR USE


The brace is indicated for combined and non-combined instabilities of the:

- Anterior cruciate ligament (ACL)
- Medial collateral ligament (MCL)
- Lateral collateral ligament (LCL)
- Posterior cruciate ligament (PCL)

### FEATURES AND BENEFITS

- Cold-malleable aluminum frame, trimmable straps, and field-serviceable
- D-rings accommodate various patient anatomies
- Step-through design moves the rigid frame cuff behind the calf
- Tool-free polycentric hinge with ROM stop options for range-of-motion restriction
- Össur Hinge Alignment Guide, Blue is You™ touch points, and tactile numbered strap tabs guide the fitting process
- Optional PCL Kit accessory for multi-ligament support

### REQUIRED MEASUREMENTS



Measure knee caliper medio-laterally at joint width in standing, weight-bearing position.



Measure thigh circumference 15 cm (6") above mid-patella. Leg in extension.

### SMARTMEASURE®

A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace. Utilizes a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.



Description	Part#
SmartMeasure Measuring Kit	CMS0000
Constant Force Caliper	B-312000000

Scans, CMS (CTI only), CCS and casting also accepted. Knee M/L required

# REBOUND® DUAL ST

## ORDERING INFORMATION


### REBOUND DUAL ST

Size	Caliper	Circumference	Right	Left
XS	3-3.5" (7.6-8.9 cm)	15-16.5" (38.1-41.9 cm)	B-242520001	B-242620001
S	3.5-4" (8.9-10.2 cm)	16.5-18.5" (41.9-47 cm)	B-242520002	B-242620002
M	4-4.5" (10.2-11.4 cm)	18.5-20.5" (47-52.1 cm)	B-242520003	B-242620003
L	4.5-5" (11.4-12.7cm)	20.5"-22.5" (52.1-57.2 cm)	B-242520004	B-242620004
XL	5-5.5" (12.7-14cm)	22.5"-24.5" (57.2-62.2 cm)	B-242520005	B-242620005
XX	5.5-6" (14-15.2cm)	24.5"-27.5" (62.2-69.9 cm)	B-242520006	B-242620006
3XL	6-6.5" (15.2 – 16.5 cm)	27.5"-30" (69.9-76 cm)	B-242520007	B-242620007

### REBOUND DUAL ST CUSTOM

Description	Right	Left
Custom	B-155520001	B-155620002

### CUSTOM BRACE COLORS



## ACCESSORIES

### OFF-THE-SHELF LINER KITS

Size	Liner Kit, Black*
XS	B-701090120
S	B-701090121
M	B-701090122
L	B-701090123
XL	B-701090124
XXL	B-701090125
3XL	B-701090126

### CUSTOM LINERS

Size	Calf, Active Grip, Black	Thigh, Doeskin, Black
XS	B-701090180	B-701090170
S	B-701090181	B-701090171
M	B-701090182	B-701090172
L	B-701090183	B-701090173
XL	B-701090184	B-701090174
XXL	B-701090185	B-701090175
3XL	B-701090186	B-701090176

\*Includes Doeskin thigh liner and Sensil calf liner

### LEGACY OFF-THE-SHELF LINERS

Size	Calf, Active Grip, Silver	Thigh, Doeskin, Silver
XS	B-701090074	B-701090067
S	B-701090075	B-701090068
M	B-701090076	B-701090069
L	B-701090077	B-701090070
XL	B-701090078	B-701090071
XXL	B-701090079	B-701090072
3XL	B-701090080	B-701090073

# REBOUND® DUAL ST

## STRAP KITS & PAD KITS

Size	Replacement Strap Kit	PCL Strap Kit
XS-L	B-701090107	B-701090065
XL-3XL	B-701090108	B-701090066

## ACTIVE GRIP STRAP PAD KIT

Size	Part#
XS-M	B-701090039
L-XXL	B-701090040

## SILICONE STRAP PAD KIT

Description	Part#
Pad Kit	B-701090037

## OTHER

Description	Part#
*Flexion/Extension Stop Set - Tool-Free Hinge	B-701090088
**Flexion Stop Set - Custom & Legacy Hinge	B-701090004
***Extension Stop Set - Custom & Legacy Hinge	B-701090008
D-Ring Kit	B-701090059
Condyle Pad Kit	B-701090009
Alligator Tab 1	RDP61106
Alligator Tab 2	RDP61102
Alligator Tab 3	RDP61103
Alligator Tab 4	RDP61104
Alligator Tab 5	RDP61105

\*Extension Stops: 0°,5°,10°,20°,30°,45° Flexion Stops: 0°,10°,20°,30°,45°,60°,75°,90°

\*\*Flexion Stops: 0°,30°,45°,60°,75°,90°

\*\*\* Extension Stops: 5°,20°,30°,40°

## SLEEVES

Size	Thigh Circumference	Sofsleeve	Neoprene Undersleeve	Neoprene Oversleeve	Padded Sports Oversleeve
S	< 15" (38 cm)	20643	20251	20051	20448
M	15 – 18" (38 – 46 cm)	20644	20252	20052	20448
L	18 – 19.5" (46 – 50 cm)	20644	20253	20053	20449
XL	19.5 – 21" (50 – 55 cm)	20645	20254	20054	20449
XXL	21.5 – 24.5" (55 – 62 cm)	20645	20255	20055	N/A
3XL	24.5 – 28" (62 – 71 cm)	20646	20256	20056	N/A
4XL	28 – 34" (71 – 86 cm)	20647	20257	20057	N/A

Size	Thigh Circumference	Sleeve with Mesh*
XS	<17" (43 cm)	SLV60011
S	17 – 18" (43 – 45 cm)	SLV60012
M	18 – 20" (45 – 51 cm)	SLV60013
L	20 – 21" (51 – 53.5 cm)	SLV60014
XL	21 – 22" (43.5 – 56 cm)	SLV60015
XXL	22 – 24" (56 – 61 cm)	SLV60016

Size	Thigh Circumference	Össur Knee Undersleeve**
XS	12 – 15" (30.5 – 38 cm)	B-250150101
S	15 – 18" (38 – 45.5 cm)	B-250150102
M	18 – 22" (45.5 – 56 cm)	B-250150103
L	22 – 26" (56 – 66 cm)	B-250150104
XL	26 – 29" (66 – 73.5 cm)	B-250150105
XXL	29 – 33" (73.5 – 84 cm)	B-250150106

\*Sleeve with Mesh: Main fabric: 92% Cotton/8% Spandex - Mesh: 100% Polyester - Band: 80% Polyester/20% Spandex

\*\*Össur Knee Undersleeve: 80% Nylon/20% Spandex