

REBOUND® HIP

Össur's Rebound Hip brace has a low-profile, lightweight design that supports patient comfort and compliance. It delivers gentle controlling forces and effective functional ROM restriction in extension/flexion and abduction/adduction to promote optimal post-surgical outcomes. The universal frame and softgood are easy to stock, assemble and quickly fit on patients.



INDICATIONS FOR USE

Hip conditions that could benefit from motion restriction following surgery:

- Femoral Acetabular Impingement arthroscopic surgery
- Gluteus medius repair surgery
- Hamstring repair surgery

FEATURES AND BENEFITS

- Low-profile, lightweight design that supports patient comfort and compliance
- Quick-release pegs easily adjust the desired range-of-motion needed for the patient's optimal rehab
- Quick-fit with sliding mechanism on frame to adjust for patient's height
- Optional, adjustable ROM hinge from -10°, 0° 10°, or 20°

REQUIRED MEASUREMENTS

- Circumference at iliac crest height
- Length from greater trochanter to mid patella

ORDERING INFORMATION

REBOUND HIP

Description	Part#
Rebound Hip	98013XRB
Rebound Hip with adjustable hinge	98012XRB
Rotation Strap	RHK-100024

X = 0 (Preset to left) or 5 (Preset to right)

*Frame: Universal; Softgood: Universal fitting to 48" hip circumference

ACCESSORIES

Description	Part#
OR Kit	980301
Belt Extender (Fitting 49-67" hip circumference)	980121RB
Fitting Tool	B-040100000