

## The Swiss Army Knife of Orthopedic Boots

VACOcast Fit is the most versatile boot on the market with **3-in-1 usability**. It's fixed plantar-flexion position and patented *Vacuum Relief* insert help patients heal faster while providing superior stability and comfort.

VACOcast's *Vacuum Relief* insert **works by sucking all of the air out** to allow microbeads to perfectly, and gently, mold to the limb. With no inward pressure, VACOcast can be worn immediately post-surgery. Our "perfect fit" can't be matched, **just try one on and experience the difference of vacuum.**

[www.vacocast.com](http://www.vacocast.com)

### 3-in-1 Usability

#### Splint

Creates a rigid structure within seconds, to secure the limb.

#### Cast

Forms a sturdy cast in less than 3-minutes with no waste.

#### Walker

Stability and style work in unison for patients to walk in comfort.



### It sucks!

That's how it works.



In this example, air is forced out of a container to form a perfect mold of the hand. VACOcast uses this same principle to create a perfect fitting orthopedic boot, unmatched in its stability and comfort by any other boot on the market.

See it in action at [www.vacocast.com](http://www.vacocast.com)